


MASTER LABEL

DATE PRINTED

31/106

PRINTED BY:

DM-8



Hồ Thị Uyên Khuong
Country Manager

AUTOPAK LABEL

**Modified Thayer –
Martin (MTM II) Agar**
221567
(4321567)
LINE# 9 0001

Approximate Formula* per Liter of Purified H₂O

Pancreatic Digest of Casein 7.5g, Selected Meat
Peptone 7.5g, Corn Starch 1.0g, Dipotassium
Phosphate 4.0g, Monopotassium Phosphate
1.0g, Sodium Chloride 5.0g, Agar 12.0g,
Hemoglobin 10.0g, Trimethoprim Lactate
5.0mg, IsoVitalex™ Enrichment 10.0mL,
V – C – N Inhibitor 10.0mL

In vitro Diagnostic

LOT 1234543 : 2004 12 12
(YYYY MM DD)



(01)30382902215678



(10)1234543(17)041212(30)0020

MADE IN USA **PRE – PROD** L – 000707(A)

LABEL ROOM LABEL

**Modified Thayer –
Martin (MTM II) Agar**
221567
(4321567)
0044

Approximate Formula* per Liter of Purified H₂O

Pancreatic Digest of Casein 7.5g, Selected Meat
Peptone 7.5g, Corn Starch 1.0g, Dipotassium
Phosphate 4.0g, Monopotassium Phosphate
1.0g, Sodium Chloride 5.0g, Agar 12.0g,
Hemoglobin 10.0g, Trimethoprim Lactate
5.0mg, IsoVitalex™ Enrichment 10.0mL,
V – C – N Inhibitor 10.0mL

in vitro Diagnostic

LOT 1234543 : 2004 12 12
(YYYY MM DD)



(01)30382902215678



(10)1234543(17)041212(30)0020

MADE IN USA L – 000707(A)

Nhãn phụ

Môi trường chuẩn bị sẵn

Mã sản phẩm: 221567

Số lưu hành:

Số lô: Xem “**LOT**” trên nhãn gốc.

Hạn sử dụng: Xem “**⌚**” trên nhãn gốc.

Hãng sản xuất/ Địa chỉ: Becton Dickinson and Company (BD), BD Diagnostic Systems/ 250 Schilling Circle Cockeysville, Maryland 21030, Mỹ.

Xuất xứ: Mỹ

Chủ sở hữu số lưu hành: Văn phòng đại diện Becton Dickinson Asia Limited tại TP. Hồ Chí Minh.

Địa chỉ: Phòng 1106-07, tầng 11, toà nhà Mê Linh Point, số 2 Ngô Đức Kế, Phường Bến Nghé, Quận 1, thành phố Hồ Chí Minh.

Nhập khẩu & phân phối bởi:

Địa chỉ:

Cảnh báo, hướng dẫn sử dụng, hướng dẫn bảo quản: Xem trong *Hướng dẫn sử dụng* đính kèm HOẶC *Website*:

TM. VPĐD Becton Dickinson Asia Limited tại TP. Hồ Chí Minh



Hồ Thị Uyên Khương

TM. Trưởng đại diện