

MASTER LABEL

DATE PRINTED 3/4/06

PRINTED BY: PM-Set



Ho Thi Uyen Khuong
Country Manager

AUTOPAK LABEL

221568
(4321568)
LINE# 9 0002
Modified Thayer -
Martin (MTM II) Agar


Approximate Formula* per Liter of Purified H₂O
Pancreatic Digest of Casein 7.5g, Selected Meat
Peptone 7.5g, Corn Starch 1.0g, Dipotassium
Phosphate 4.0g, Monopotassium Phosphate
1.0g, Sodium Chloride 5.0g, Agar 12.0g,
Hemoglobin 10.0g, Trimethoprim Lactate
5.0mg, IsoVitalex™ Enrichment 10.0mL,
V-C-N Inhibitor 10.0mL

in vitro Diagnostic

LOT

2222222

2004 12 12
(YYYY MM DD)



(01)30382902215685



(10)2222222(17)041212(30)0100

POST-PROD L-000708(A)

MADE IN USA

LABEL ROOM LABEL

221568
(4321568)
0045
Modified Thayer -
Martin (MTM II) Agar

Approximate Formula* per Liter of Purified H₂O
Pancreatic Digest of Casein 7.5g, Selected Meat
Peptone 7.5g, Corn Starch 1.0g, Dipotassium
Phosphate 4.0g, Monopotassium Phosphate
1.0g, Sodium Chloride 5.0g, Agar 12.0g,
Hemoglobin 10.0g, Trimethoprim Lactate
5.0mg, IsoVitalex™ Enrichment 10.0mL,
V-C-N Inhibitor 10.0mL

in vitro Diagnostic

LOT 2222222

2004 12 12
(YYYY MM DD)



(01)30382902215685



(10)2222222(17)041212(30)0100

L-000708(A)

MADE IN USA

Nhãn phụ

Môi trường chuẩn bị sẵn

Mã sản phẩm: 221568

Số lưu hành:

Số lô: Xem “**LOT**” trên nhãn gốc.

Hạn sử dụng: Xem “**⌚**” trên nhãn gốc.

Hãng sản xuất/ Địa chỉ: Becton Dickinson and Company (BD), BD Diagnostic Systems/ 250 Schilling Circle Cockeysville, Maryland 21030, Mỹ.

Xuất xứ: Mỹ

Chủ sở hữu số lưu hành: Văn phòng đại diện Becton Dickinson Asia Limited tại TP. Hồ Chí Minh.

Địa chỉ: Phòng 1106-07, tầng 11, toà nhà Mê Linh Point, số 2 Ngô Đức Kế, Phường Bến Nghé, Quận 1, thành phố Hồ Chí Minh.

Nhập khẩu & phân phối bởi:

Địa chỉ:

Cảnh báo, hướng dẫn sử dụng, hướng dẫn bảo quản: Xem trong *Hướng dẫn sử dụng* đính kèm HOẶC *Website*:

TM. VPĐD Becton Dickinson Asia Limited tại TP. Hồ Chí Minh



Hồ Thị Uyên Khương

TM. Trưởng đại diện