




VA-8002
Neonate ~ Adult

REUSABLE

Test Lung w/1L Silicone Bag

1. Tidal Volume max 1000ml
2. Autoclavable temperature 120 degree C
3. Material Silicone
4. Standard Resistance:
Rp 50cmH2O/L/S for TV 0ml~200ml



LOT 607231

2020-08



Great Group Medical Co., Ltd.
No 168, Xingong 2nd Rd, Changhua County, Taiwan

ECREP Wellkang

Suite B, 29 Harley Street, London W1G9QR England, UK