

Test Report

Name of the sample: x-ray shielding lead glass

Number of the sample: 2016-327

Dimension:1000mmx800mmx10mm ZF3

Quantity of the sample: 1 piece

Description of sample: light yellow, diaphanous solid

Examination item: shielding capability of x-ray protective material

Examination standards: attenuation property test of x-ray protective material (GBZ/t147-2002)

Main equipment and the number: NE2550 the second degree doser J-067

**DCI8500 precise current integrator TK30
ionization chamber J-102**

Date of examination:2016.04.04

date of sending the sample :2016.04.01

Name of applicant company:

JiaoZuo Changxin Technical Development of Radiation Protection co.,Ltd

Address of Applicant Company: No.3 Jiankang Road

Zip Code of Applicant Company: 454150

Tel: 0391-2927111

Report of the test:

Lead equivalency: 2.04mmpb (120kv 2.5mmAL)

Expanded uncertainty: 6.0 % (k=3)

检测报告

第 1 页, 共 1 页

样品名称: 铅玻璃板

样品编号: 2016-327

型号规格: ZF3 1000mm×800mm×10mm

样品数量: 一块

样品性状: 淡黄色透明固体板状

检测项目: X 射线防护材料屏蔽性能检测

检测依据: X 射线防护材料衰减性能的测定 (GBZ/T147-2002)

主要仪器设备及其编号: NE2550 二级标准剂量仪 J-067

DC18500 精密电流积分仪 TK30 电离室 J-102

检测日期: 2016 年 4 月 4 日 送样日期: 2016 年 4 月 1 日

委托单位名称: 焦作市长信防辐射技术开发有限公司

委托单位地址: 焦作市健康路 3 号

委托单位邮编: 454000

联系电话: 15139168880


检测结果:

铅当量: 2.4mmPb (120 kV 2.5 mA1)

扩展不确定度: 6.0% ($k=3$)

以下空白

授权签字人:



签发日期: 2016年4月4日

