

Crown Removers

As instrument to remove crown hard to remove easily, first, make gap in crown using bur and twist after putting blade of instrument in gap.

Prosthodontic
Crown Removers

Crown Remover

NEW PRODUCT
3RECDA

· For removing troublesome crowns with relative ease. The 5 etching points make it possible to access to a crown from various sides.

 Autoclavable



NEW PRODUCT
3RECDB

· Same instruction as 3RECDA but has one point

 Autoclavable



RECR3KIT

· RECR3-A : To put and pull gap between Pontic and abutment of bridge.
· RECR3-B : To hang and pull sawtooth of bottom of margin of single crown
· RECR3-C : To hang and pull hook of temporary single crown.

RECR3-A

Crown Remover Tip
· Tip A



RECR3-B

Crown Remover Tip
· Tip B



RECR3-C

Crown Remover Tip
· Tip C



PRODUCTS FOR DENTISTRY
OSUNG MIND CO.,LTD.

Crown Forceps · Crown Gripper


Crown Forcep

· Easy to remove crown with rubber on both tips.

NEW PRODUCT

CF01

· Length : 145mm(±5mm)
· Rubber Tip 20pcs, Powder 5g



NEW PRODUCT

CF02

· Length : 145mm(±5mm)
· Rubber Tip 20pcs, Powder 5g



NEW PRODUCT

CFRT85

· 20pcs
· Rubber Tip



NEW PRODUCT

CFP5

· 5g
· Powder



Crown Gripper

· Easy to remove crown & temporary tooth.

NEW PRODUCT

CG01

· Length : 155mm(±5mm)



Prosthodontic
Crown Forceps ·
Crown Gripper

PRODUCTS FOR DENTISTRY
OSUNG MIND CO.,LTD.