

Silicone Wound Drainage System



Indication

Reservoirs are intended to be used for post-surgery for collecting the exudate.

Feature

- Anti-back valve to prevent bacteria from reflex of fluid.
- An easy-clip is provided to attach to bed, wheelchair, or IV stand for mobility.
- 100% medical grade silicone meets biocompatibility requirement.
- Low vacuum design provides continue gentle suction and minimizes tissue trauma.
- Size available: 100ml, 150ml, 200ml, 400ml, top-outlet, and bottom-outlet.
- Bottom-outlet is designed to connect with collecting bag.
- Sterile double package.



Order Information

Code	Description
CW-69101	Top Outlet 100cc
CW-69201	Top Outlet 150cc
CW-69303	Top Outlet 200cc
CW-69402	Top Outlet 400cc (Rectangular)
CW-69404	Top Outlet 400cc
CW-69224	Bottom Outlet 150cc with Collecting Bag
CW-69353	Bottom Outlet 200cc with Collecting Bag



400cc (Rectangular)



400cc

Collecting Bag

Feature

- Various types of connector for easy application.
- PVC 500ml bag, easy for connection.
- Sterile double package.

Order Information

Code	Description
CW-69801	2.4 mm connector
CW-69802	5.4 mm connector
CW-69803	Adjustable connector
CW-69804	Female Universal connector
CW-69805	4 mm connector

Connector Type



2.4mm Connector

5.4mm Connector

Adjustable Connector

Female Universal Connector

4mm Connector

Connector



Drainagebeutel 500ml
Patient name
Date
Collection/Drainage

500ml
400
300
200
100
50