



# Heating Pads

## Part Number

- BE-220 Heating Pads

## Features

- Soothing warmth relieves tired, aching muscles and temporary arthritic pain.
- Besmed heating pads are versatile and light weight, comfortable to use.
- Use them on most parts of the body, arms legs, shoulders, back abdomen, elbows, knees, ankles.
- Absorbent sponge provides penetrating moist heat.
- Soft cover with tie bands to hold the heating pad in place.
- Removable, washable cover with snap closure.
- 100% moisture-proof, vinyl covered pad.

## Specifications

- Heat type: dry/moist
- Size: 30x39cm
- VOLTAGE: 110V & 220V

## Packaging

- BE-220 12 Pcs/Case
- N.W. / CTN: 17.56 KG
- G.W. / CTN: 18.47 KG
- CBM / CTN: 0.05
- PACKAGE: BOX