



No: 8421/PKN-VKNQG

TEST REPORT

1. Name of sample: 4-layer antimicrobial medical mask, Brand: Phuchthanhmedical.
Manufactured at Kim Bang Pharmaceutical Joint Stock Company.
Address: No. 20 Quang Trung Street, Que Town, Kim Bang District,
Ha Nam
2. Sample code : 04202536/DV.1
3. Sample description: Sample is put in a paper box, 50 pieces/box - Quantity: 1
MFG - EXP: No information; No storage sample
4. Number of sample: 01 sample
5. Sample storage time: No storage
6. Sample received date: 21/04/2020
7. Testing time: 21/04/2020 - 29/04/2020
8. Customer: Kim Bang Pharmaceutical Joint Stock Company
Address: No. 20 Quang Trung Street, Que Town, Kim Bang District,
Ha Nam
9. Test result(s): Microbiological test

| No. | Parameter | Unit | Test method | Result | Compare with BS EN 14683:2019 |
|-----|---------------------------------|--------|---|---------|-------------------------------------|
| 9.1 | Bacterial filtration efficiency | /piece | Ref. TCVN 8389-2:2010 & BS EN 14683:2019 | 99.99 % | Satisfactory |
| 9.2 | Microbial cleanliness | CFU/g | BS EN 14683:2019 | 13.75 | Satisfactory |

Note: Bacterial filtration efficiency is equivalent to the protective ability against bacterial infection of the face mask

Hanoi, April 29th, 2020

On behalf of Director General
Head of Quality Assurance Laboratory



Le Thi Phuong Thao, PhD.