

The fully-functional, streamlined arthroscopy cart



Manufactured by:



Distributed by:



Big features, small footprint

This sleek and simply-designed arthroscopy cart provides full functionality and many unique features without the space required by larger carts. And it is customizable with a three or four shelf configuration to suit the specific needs and requirements of your OR.

Sleek design

with minimal gaps helps prevent germ formulation and makes cleaning easy

Lightweight and durable

construction can accommodate up to 176lbs/80kg



Hollow bearing column facilitates cord compartmentalization and containment

On/off switch can be installed at top of support column

Ergonomic handle can be set to 90° for easy moving or 15° to minimize footprint



Fold-out cover means easy access to data lines and power cables

Customizable construct features three or four shelf configurations

Rubber lips permit cable entry and exit over the entire height of the column



Convenient bottom drawer offers additional storage

Wide base provides stability and facilitates easy handling

Ordering Information

Arthroscopy Cart

Reference #	Description
72205119	3 Shelf Arthroscopy Cart
72205120	4 Shelf Arthroscopy Cart
72205121	Replacement Shelf
72205124	CO ² Tank Holder

Manufactured by:
ITD GmbH

Distributed by:

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810 USA

www.smith-nephew.com
T +1 978 749 1000
US Customer Service: +1 800 343 5717

Isolation transformer protects staff and equipment from electrical surge