







Bring safety & belief to your test!

QUE CẤY NHỰA

STICKS/LOOPS





QUE CẤY NHỰA/LOOPS

Que cấy nhựa Midamec được sản xuất bằng nhựa HIPS chánh phẩm, chất liệu mềm dẻo dễ dàng lấy mẫu cấy vi sinh mà không ảnh hưởng đến môi trường cấy.

Đầu que cấy vi sinh tròn, không bavia, lý tưởng cho các quy trình cấy vi sinh, lấy mẫu đạt độ chính xác cao.

Đóng gói kín, không thấm nước ngăn chặn hơi ẩm ô nhiễm.

Vô trùng bằng tia gamma.

Được sản xuất, đóng gói theo tiêu chuẩn phòng sạch cấp độ ISO 5, theo tiêu chuẩn ISO 9001: 2008, ISO 13485: 2003

Plastic loop is perfect circle & made by flexible virgin HIPS that does not cause to the culture environment.

MidaMec Loops are calibrated for culturing, transferring, streaking or isolating tissue cultures in cytology or histology.

Package in PP zip bag.

Sterile by gamma ray, especially perfect pure.

Under Self-contained Manufacturing and packing in clean room environment class ISO5, complied with ISO 9001:2008, ISO 13485: 2003

THÔNG SỐ KỸ THUẬT / SPECIFICATIONS

Mã sản phẩm/ Product Code	Nhựa/ Resin	Sản phẩm/ Product Name	Dung tích/ Volume	Màu/ Color	Đóng gói/ Packing	Vô trùng/ Sterilization
LP0001-2	HIPS	Que cấy 1ul/ Stick/Loop 1ul	1ul	Xanh/ Blue	2 pcs/pouch 8,000 pcs/carton	Theo yêu cầu bằng tia gamma/ As requirement by Gamma
LP0001-10	HIPS	Que cấy 1ul/ Stick/Loop 1ul	1ul	Xanh/ Blue	10 pcs/bag 8,000 pcs/carton	
LP0001-25	HIPS	Que cấy 1ul/ Stick/Loop 1ul	1ul	Xanh/ Blue	25 pcs/bag 8,000 pcs/carton	
LP0001-100	HIPS	Que cấy 1ul/ Stick/Loop 1ul	1ul	Xanh/ Blue	100 pcs/bag 8,000 pcs/carton	
LP0002-2	HIPS	Que cấy 10ul/ Stick/Loop 10ul	10ul	Vàng/ Yellow	2 pcs/pouch 8,000 pcs/carton	
LP0002-10	HIPS	Que cấy 10ul/ Stick/Loop 10ul	10ul	Vàng/ Yellow	10 pcs/bag 8,000 pcs/carton	
LP0002-25	HIPS	Que cấy 10ul/ Stick/Loop 10ul	10ul	Vàng/ Yellow	25 pcs/bag 8,000 pcs/carton	
LP0002-100	HIPS	Que cấy 10ul/ Stick/Loop 10ul	10ul	Vàng/ Yellow	100 pcs/bag 8,000 pcs/carton	