



Bring safety & belief to your test!

TUBE VI SINH & ĐÔNG LẠNH

MICRO TEST TUBES & CRYOTUBES





TUBE VI SINH & ĐÔNG LẠNH/MICRO TEST TUBES & CRYOTUBES

Được dùng để lưu trữ trong tủ đông.

Quay ly tâm tách mẫu.

Được sản xuất bằng nhựa PP chánh phẩm.

Thành ống dày chịu được quay ly tâm 18,000 RPM

Sản phẩm trong, dễ dàng quan sát.

Được thiết kế liền thân, đạt độ kín cao khi đóng nắp.

Có thể viết lên ống hoặc nắp.

Sản phẩm có thể đun trong nước sôi.

Use for cold store.

Use in centrifugation extract specimen.

Made from virgin PP plastic.

Thick wall for centrifugation up to 18,000 RPM.

Extra clarity for better visual inspection.

One piece designed - Snug-fitting cap.

Frosted writable on closure and sides of tubes.

Boil-proof design.

THÔNG SỐ KỸ THUẬT / SPECIFICATIONS

Mã sản phẩm/ Product Code	Nhựa/ Resin	Sản phẩm/ Product	Dung tích/ Volume	Màu/ Color	Đóng gói/ Packing	Vô trùng/ Sterilization
CT0001	PP	Tuýp vi sinh 1.5ml/ Eppendorf 1.5ml	1.5ml	Tự nhiên/ Natural	1,000pcs/PE bag 7,000pcs/carton	Theo yêu cầu/ As requirement