



Beckman Coulter Ireland Inc.,
Lismeehan,
O'Callaghan's Mills,
Co. Clare,
Ireland.
Tel. +353-65-6831100.
Fax. +353-65-6831122.

Technical Specification

iChem VELOCITY Wash Solution

Part Number: 800-7217

- Intended use:** The iChemVELOCITY Wash Solution is a buffer solution for routine rinse of iChemVELOCITY Urine Chemistry System. This reagent is intended for in vitro diagnostic use with the iChemVELOCITY.
- Packaging**
 - 2 x 7L iChem VELOCITY Wash Solution
 - 2 x 1 μ m Filter
- Storage:** 20 – 28°C
- Physical and Chemical Specification**

Attribute	Specification
Physical State	Liquid
Color	Colorless
Odor	Aromatic odor
OD (450 nm)	≤ 0.0100
OD (550 nm)	≤ 0.0070
OD (650 nm)	≤ 0.0070
Clarity	≤ 10 NTU
SG	1.000 ± 0.002

Authorized Signatory & Date:

Marguerita Sweeney 27th June 2019
Marguerita Sweeney
Regulatory Affairs Manager
Beckman Coulter Ireland Inc.
Phone: 353 65 6831495

Registered Office: Lismeehan,
O'Callaghan's Mills,
Co. Clare.

Directors: Piers Devereux
Paul Henry
Orlaith A. Lawler

Incorporated in Republic of Ireland
with limited liability.
Co. Reg. No. 501471.