

Công ty TNHH HOYA LENS VIỆT NAM - CHI NHÁNH QUẢNG NGÃI  
 Sản phẩm Xuất Khẩu và Không bán tại thị trường Việt Nam




Mã sản phẩm MẪU NHÃN TRÒNG KÍNH XUẤT KHẨU (THAM KHẢO)

HO150SB3 BHV65  
 HO150SB3 BHV70

**Finish Lens Label Format**

**Definition:**

- no.1 OPC Code Barcode
- no.2 2D Barcode
- no.3 Trace Abilities Barcode
- no.4 Power of Lens
- no.5 Label Name
- no.6 Lot Number
- no.7 Worker Number  
 [=100000-Actual worker no.]  
 [Ex. 100000-79973=20027]
- no.8 Packing machine number
- no.9 Packing Date [DD MM Y]  
 Ex1. 09709 => 9 July 2009  
 DD= 09, MM=70, Y=9  
 Ex2. 62219 => 26 December 2009
- no.10 Factory [HOLV=4]
- no.11 Packing Sequence number
- no.12 Diameter
- no.13 Actual CT
- no.14 Install Barcode  
 =4991696 + Turn Over + check digit  
 \* In case of JAN code
- no.15 Jan Code Barcode
- no.16 Quarter (Time period)



**Packing Quarter**

Year	Quarter
Year 2 digits + Quarter 1 digit	
W = 2005	1 = January - March
X = 2006	2 = April - June
Y = 2007	3 = July - September
Z = 2008	4 = October - December

WW = 2009  
 XX = 2010  
 YY = 2011  
 ZZ = 2012  
 W = 2013  
 Y = 2015  
 Z = 2016  
 WW = 2017

1.50

HL15BR3S-H

SINGLE VISION  
 PLASTIC LENS  
 Refractive Index: 1.50

**HILUX**



**HV** HI-VISION

**Brown**  
**SUNTECH**

**HOYA**