



ISOLATION GOWN

Disposable medical protective clothing

- *For single use only
- *Store in a cool, dry place
- *Away from heat and direct sunlight



Specifications

| | |
|------------|----------------------------------|
| Product | Non-woven isolation gown and cap |
| Size | Universal, one size fits most |
| Dimensions | Length 45", width 55" |
| Color | White |
| Wrist | Knit Cuf |
| Material | Non woven |
| Pack | 2 set/bag |

MAI LAM FASHION COMPANY LIMITED
ADD: 103 Nguyen Dinh Chieu, Ward 6, District 3, HCMC
Phone: (84 8) 3520 0066