

● YFC361L Three Function Manual Bed



Size: 2160*960*500-750mm (L*W*H)

Functions:

1. HI-LOW: 500-750mm
2. Back section adjustable from 0° to 75°
3. Knee section adjustable from 0° to 40°

Structure:

Electrophoresis and powder coating
Wholly molded cold steel sleeping platform

Accessories:

Head and Foot board: removable, ABS material, detachable
Side rails: Two aluminum alloy side rails, collapsible
Castor: four 125mm castor and two brakes
I.V. pole
Maximum loading capacity: 250kgs

■Note:Color is optional

● Bed Mattress



Size according to the sleeping surface.
Full sponge mattress with PU cover , waterproof
Thickness: 15cm .
The sponge density : 35kg/m3
Four section, each with a zipper, removable and easy cleaning.

