



E003:

Refrigerators and freezers

PQS code:	E003/017
Type of appliance:	Two-mode vaccine refrigerator OR waterpacks freezer
Manufacturer's reference:	TCW 3000 AC
Manufactured in:	Luxembourg
Company:	B Medical Systems Sarl (Dometic)
Address:	17 Op der Hei L-9809 Hosingen Luxembourg
Telephone:	+ 35 2 92 07 311
Email:	info@bmedicalsistemas.com
Web address:	www.bmedicalsistemas.com

Specifications

Climate zone:	Hot	Min rated ambient temp:	+15°C
Refrigerant:	R134a	Energy source:	Electric Mains 220-240V 50/60Hz or Electric Mains 115/60Hz optional
Appliance tested at:	+43°C	Ext dimensions (HxLxD)	91 x 127 x 78 cm
Performance at:	+43°C	Fuel and cycle type:	Electric - compression

REFRIGERATOR

FREEZER

Freeze protection: Grade A

Vaccine storage capacity: (Litres)	150		
Gross volume (Litres):	204	Gross volume:	204
Holdover time (hours):	53 hours 10 min; 3 hours 55 min	Waterpack freezing capacity:	N/A
Energy consumption, stable running (kWh/24 hours):	1.37 kWh/24h; 3.83 kWh/24h	Waterpack storage capacity:	112.2 kg
Energy consumption, cool down test (kWh/24 hours):	1.32 kWh/24h; 3.91 kWh/24h	Energy consumption during freezing:	-----