

Tên hàng : **Dụng cụ kiểm tra huyết áp (người lớn)**

Quy cách : **hộp/cái**

Xuất xứ : **Tanaka - Nhật**

Lô SX và HD: **không hạn dùng**

NNK & PP: **Cty CP Trang Thiết Bị Kỹ Thuật Y Tế TPHCM**

- 252-254 Cống Quỳnh, P. Phạm Ngũ Lão, Q.1, TPHCM

- ĐT : 02838325712



<Specifications>

- *Manometer : 20 – 300mmHg
- *Accuracy : ± 3 mmHg
- *Air System : Velcro Cuff with Latex Bag
: Bulb with Valves
- *Weight : 430g
- *Packing : 50sets/carton
: 0.065m³, 23.0kg

Tên hàng : **Dụng cụ kiểm tra huyết áp (trẻ em)**

Quy cách : **hộp/cái**

Xuất xứ : **Tanaka - Nhật**

Lô SX và HD: **không hạn dùng**

NNK & PP: **Cty CP Trang Thiết Bị Kỹ Thuật Y Tế TPHCM**

- 252-254 Cống Quỳnh, P. Phạm Ngũ Lão, Q.1, TPHCM

- ĐT : 02838325712



1

Cert Certification is the notation of No.500-C

1 **Aneroid Sphygmomanometer for Pediatrics No.500-C1 ~ C5**

- *Manometer : 20 ~ 300mmHg
(No.500-M) : with zero stop pin
- *Accuracy : ± 3 mmHg
- *Air System : Velcro Cuff for Pediatrics
with Latex Bag
: Bulb & Valves (No.90)

Item No.	Cuff Size	Type
No.500-C1	4.0 × 26.0 cm	New Born
No.500-C2	4.5 × 26.0 cm	Infant
No.500-C3	6.5 × 28.0 cm	Children
No.500-C4	8.5 × 30.0 cm	Children-M
No.500-C5	10.5 × 34.0 cm	Children-L

- *Weight : 300 g(No.500-C3)
- *Packing : 50 sets/carton
: 0.065m³, 16.0 kg