





## 1. M-30R RINSE

<b>mindray</b>		<b>M-30R RINSE</b>		<b>5.5L</b>	
				P/N:A12-000197- - -	
<p><b>For use on MINDRAY BC-3000Plus/BC-3000CT/BC-3200/BC-3300/BC-3600 hematology analyzer. For use as an azide-free, filtered isotonic solution for cleaning hematology analyzer. Use as directed in instrument operation manual.</b></p>				<p><b>IVD</b> For In Vitro Diagnostic Use</p>	
<p><b>WARNINGS:</b> Do not inhale and/or ingest. In case of inhaling or ingestion, immediately go see a doctor. Avoid skin contact. In case of skin contact, flush the affected area with plenty of water. Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and go see a doctor.</p>		<p><b>ACTIVE INGREDIENTS</b> Sodium Sulfate Anhydrous..... 3.0-5.5 g/L Sodium Chloride ..... 7.5-11.5 g/L Anti-fungal and Anti-bacterial Agents .... 1.0-3.0 g/L Polyoxyethylene Ester ..... 0.3-1.5 g/L</p>		<p>30°C 2°C</p> <p><b>CE</b></p>	
<p><b>STORAGE AND HANDLING:</b> Store product at 2°C to 30°C(35°F to 86°F). Keep container closed. If product has been frozen(either partially or completely), let product warm to room temperature before mixing. Mix product by gentle inversion. Verify background results before analyzing patient samples. If the operating environment is 15°C to 30°C(59°F to 86°F), discard opened container after 60 days. Dispose of waste product, unused product and contaminated packaging in accordance with governmental requirements.</p>		<p> <b>LOT</b></p>			
<p> SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China</p>		<p><b>EC REP</b> Shanghai International Holding Corp. GmbH(Europe) Eiffestraße 80, 20537 Hamburg Germany</p>			

## 2. M-30E E-Z CLEANSER

<b>mindray</b>		<b>M-30E E-Z CLEANSER</b>		<b>100mL</b>	
				P/N:A12-000045- - -	
<p><b>STORAGE AND HANDLING:</b> Store product at 2°C to 30°C(35°F to 86°F). If the operating environment is 15°C to 30°C (59°F to 86°F), discard opened container after 60 days. Keep container closed when not in use. Dispose of waste product, unused product and contaminated packaging in accordance with governmental requirements.</p>				<p><b>For use on MINDRAY BC-3000Plus, BC-3000CT, BC-3200 hematology analyzer. For use as a cleaning reagent for cleaning hematology analyzer periodically. Use as directed in instrument operation manual.</b></p>	
<p><b>CE</b> 2°C 30°C</p> <p><b>IVD</b> For In Vitro Diagnostic Use</p>		<p><b>WARNINGS:</b> Do not inhale and/or ingest. In case of inhaling or ingestion, immediately go see a doctor. Avoid skin contact. In case of skin contact, flush the affected area with plenty of water. Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and go see a doctor.</p>		<p><b>ACTIVE INGREDIENTS</b> Proteolytic Enzyme ..... 3.0-10 g/L Surfactant ..... 0.3-1.5 g/L Sodium Chloride ..... 3.0-5.0 g/L Anti-fungal and Anti-bacterial Agents • 0.5-2.5 g/L Buffering Agents ..... 1.0-4.0 g/L</p> <p> <b>LOT</b></p>	
<p><b>EC REP</b> Shanghai International Holding Corp. GmbH(Europe) Eiffestraße 80, 20537 Hamburg Germany</p>		<p> SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China</p>			