

24<sup>th</sup> Dec 2021

**Product:** Durex Play Warming Lube

**Formulation code:** 8032226

**Manufacturing site:** RB Healthcare Manufacturing, Thailand

Dear Sir or Madam,

We are providing you with information of Durex Play Warming Lube as following;

**Formulation specification:**

Test	Limits	
	Release Specification	Shelf Life Specification
Appearance	Clear, colourless gel free from lumps and extraneous matter	Clear, colourless gel free from lumps and extraneous matter
Odour	Slightly sweet	Slightly sweet
pH at 25 °C	3.8-4.8	3.8-4.8
Viscosity - RV4, 50 rpm at 25°C (cP)	1,500 - 4,000	Record Only
Benzoic Acid (%w/w)	0.143 - 0.157	0.130-0.174
TAMC Frequency for testing; Release spec: One batch per year Shelf life spec: As per Analysis Schedule	NMT 200 cfu/g	NMT 200 cfu/g
TYMC Frequency for testing; Release spec: One batch per year Shelf life spec: As per Analysis Schedule	NMT 20 cfu/g	NMT 20 cfu/g
Content of microorganisms; Frequency for testing; Release spec: One batch per year Shelf life spec: As per Analysis Schedule Pseudomonas Aeruginosa Staphylococcus Aureus Candida albicans Escherichia coli Preservative Efficacy Testing: (Shelf-life specification only)	Absent/g Absent/25g Absent/g Absent/g -	Absent/g Absent/25g Absent/g Absent/g Comply with EP

Justification reference no: JU138/21 v1.0

### AUTHOR AND APPROVAL

Role	Print Name	Signature	Date
Author-R&D Associate Operations Intimate Wellness	Kingfa Yaowarat	Kingfa Y.	24 - Dec - 2021
Approver- Senior R&D Associate Operations Intimate Wellness	Monthira Reodacha	MR	24 Dec 21

### REVISION HISTORY

Version Number	Date	Reason for Change	Author
1.0	24 <sup>th</sup> Dec 2021	New justification	Kingfa Yaowarat