

GIÁM ĐỐC



Trần Ngọc Mỹ

Nhãn phụ:

Thiết bị điều trị da



Chủng loại: HOLLYWOOD SPECTRA

Số lưu hành: XXXXXX

Chủ sở hữu số lưu hành: Công ty TNHH Thiết bị y tế Inamed, địa chỉ: tầng 8, phòng 8.1, tòa nhà Etown, số 364 Cộng Hòa, Phường 13, Quận Tân Bình, TP.HCM


Cơ sở bảo hành: Công ty TNHH Thiết bị y tế Inamed, địa chỉ: tầng 8, phòng 8.1, tòa nhà Etown, số 364 Cộng Hòa, Phường 13, Quận Tân Bình, TP.HCM




Chủ sở hữu/Sản xuất bởi: Lutronic Corporation, địa chỉ: Lutronic Center, 219, Sowon-ro, Deogyang-gu, Goyang-si, Gyeonggi-do, Hàn Quốc


Số seri, ngày sản xuất: xem  ,  trên nhãn gốc phía sau thân máy


Thông tin cảnh báo, hướng dẫn sử dụng, hướng dẫn bảo quản và các thông tin khác: xin xem trong Hướng dẫn sử dụng đi kèm máy


Nhãn gốc:

LUTRONIC 

Product Name : Frequency Doubled Q-switched Nd-YAG Laser
Model Name : HOLLYWOOD SPECTRA, Weight : 84.7kg
Laser Product Class : CLASS 4 
Laser Wavelength : 1064/532/585/650nm
Electrical Shock Class : CLASS I, TYPE B 
Electrical Rating : AC220-230V: 50/60Hz, 1.7kVA
AC100-120V: 50/60Hz, 1.7kVA 

 **SN**

 **Manufacturer :** Lutronic Corporation.
Lutronic Center, 219, Sowon-ro, Deogyang-gu,
Goyang-si, Gyeonggi-do, Korea
TEL : +82-31-908-3440 **FAX:** +82-31-907-3440
http://www.lutronic.com **MADE IN KOREA**

 **Before operating the device, read the operator's manual.**
EC Representative: Emergo Europe
Prinsessegracht 20, 2514 AP The Hague, The Netherlands