



宁波星创医疗器械有限公司

NINGBO SINTRUE MEDICAL INSTRUMENTS CO., LTD.

---

## Instruction for Use

**Product name:** Aneroid Sphygmomanometer Accessories

**Model:** Bladder, Bulb latex bulb with end valve

### Brief Introduction

This is Accessories for Aneroid Sphygmomanometer.

### Intended use

Used in conjunction with a Aneroid Sphygmomanometer to measure human blood pressure

### How to use:

Put the bladder into the cotton cuff and connect to the bulb with end valve.

### Notice :

Stored in the good ventilated rooms of which its temperature is lower then 30 the degree of Celsius.