



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Transparent dressing 20x10cm, with absorbent pad – Sterile**

LOT: **721746/1/A**

ARTICLE: **4381-505**

SPECIFICATION	STANDARD	RESULTS
Composition	Polyester, Polyurethane, Viscose	<i>CONFORM</i>
Backing Length	200mm ± 5mm	199mm
Backing Width	100mm ± 1mm	100mm
Woundpad Length	150mm ± 5mm	150mm
Woundpad Width	50mm ± 2mm	50mm
Pouch Length	260mm ± 7mm	262mm
Pouch Width	135mm ± 3mm	135mm
CFU/g before Steril	<10 ²	<i>CONFORM</i>

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Transparent dressing 7,2x5cm, without absorbent pad – Sterile**

LOT: **714065/1/A**

ARTICLE: **4381-510**

SPECIFICATION	STANDARD	RESULTS
Composition	Polyurethane	<i>CONFORM</i>
Backing Length	72mm ± 2,2mm	72mm
Backing Width	50mm ± 1,5mm	50mm
Pouch Length	68mm ± 3mm	69mm
Pouch Width	101mm ± 5mm	102mm
CFU/g before Steril	<10 ²	<i>CONFORM</i>

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Transparent dressing 12,5x10cm, without absorbent pad – Sterile**

LOT: **722958/1/A**

ARTICLE: **4381-511**

SPECIFICATION	STANDARD	RESULTS
Composition	Polyurethane	CONFORM
Backing Length	125mm ± 5mm	126mm
Backing Width	100mm ± 1mm	100mm
Pouch Length	185mm ± 7mm	187mm
Pouch Width	135mm ± 3mm	135mm
CFU/g before Steril	<10 ²	CONFORM

Approval:

Bastos Viegas, s.a.
Quality Management system
Sistema de gestão da qualidade
22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Orthopaedic padding, natural, 80g/m², 7,5cmx2,7m - Non Sterile**

LOT: **3B318**

ARTICLE: **420-102**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Viscose Fibres	CONFORM
Colour	White	CONFORM
Fluorescence	Absence	CONFORM
Length	2,7m ± 10%	2,67m
Width	7,5cm ± 10%	7,5cm
Weigth/ pc	16,2g ± 15%	16,4g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Orthopaedic padding, natural, 80g/m², 15cmx2,7m - Non Sterile**

LOT: **2J129**

ARTICLE: **420-104**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Viscose Fibres	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	2,7m ± 10%	2,69m
Width	15cm ± 10%	15,2cm
Weigth/ pc	32,4g ± 15%	32,8g

Approval:

[Signature]
Quality Management system
Sistema de gestão da qualidade
22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Orthopaedic padding, natural, 80g/m², 10cmx2,7m - Non Sterile**

LOT: **2K202**

ARTICLE: **420-103**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Viscose Fibres	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	2,7m ± 10%	2,7m
Width	10cm ± 10%	9,9cm
Weigth/ pc	21,6 ± 15%	21,8g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Stockinette, 5cmx25m, 100% cotton – Non Sterile**

LOT: **696527/1/A**

ARTICLE: **4213-002**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Cotton	CONFORM
Colour	Beige	CONFORM
Fluorescence	Absence	CONFORM
Length	25m ± 10%	25,2m
Width	5cm ± 10%	4,9cm
Weigth/ pc	>465,9g	468,0g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Stockinette, 7cmx25m, 100% cotton – Non Sterile**

LOT: **720137/1/A**

ARTICLE: **4213-003**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Cotton	<i>CONFORM</i>
Colour	Beige	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,4m
Width	7cm ± 10%	7,3cm
Weigth/ pc	>595,0g	598,0g

Approval:





CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Stockinette, 10cmx25m, 100% cotton – Non Sterile**

LOT: **735556/1/A**

ARTICLE: **4213-004**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Cotton	<i>CONFORM</i>
Colour	Beige	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,4m
Width	10cm ± 10%	10,2cm
Weigth/ pc	>895,4g	896,0g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Stockinette, 15cmx25m, 100% cotton – Non Sterile**

LOT: **728861/1/A**

ARTICLE: **4213-006**

SPECIFICATION	STANDARD	RESULTS
Composition	100% Cotton	CONFORM
Colour	Beige	CONFORM
Fluorescence	Absence	CONFORM
Length	25m ± 10%	25,4m
Width	15cm ± 10%	15,2cm
Weigth/ pc	>1707,9g	1709,0g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Elastic net bandage. Roll with 25m, stretched. Size: 1 – Non Sterile**

LOT: **740520/1/A**

ARTICLE: **421-001**

SPECIFICATION	STANDARD	RESULTS
Composition	70% Polyamide 30% Polyurethane	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,6m
Width	1,6cm ± 10%	1,6cm
Weigth/ pc	>45,4g	45,6g

Approval:



22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Elastic net bandage. Roll with 25m, stretched. Size: 2 – Non Sterile**

LOT: **732252/1/A**

ARTICLE: **421-002**

SPECIFICATION	STANDARD	RESULTS
Composition	70% Polyamide 30% Polyurethane	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,4m
Width	1,8cm ± 10%	1,8cm
Weigth/ pc	>56,4g	56,6g

Approval:

Alfredo
Quality Management system
Sistema de gestão da qualidade
22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Elastic net bandage. Roll with 25m, stretched. Size: 4 – Non Sterile**

LOT: **742590/1/A**

ARTICLE: **421-004**

SPECIFICATION	STANDARD	RESULTS
Composition	70% Polyamide 30% Polyurethane	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,2m
Width	2,4cm ± 10%	2,4cm
Weigth/ pc	>80,2g	80,6g

Approval:

Bastos Viegas,s.a.
Quality Management system
Sistema de gestão da qualidade
22-08-2013



CERTIFICATE OF ANALYSIS

ARTICE DESCRIPTION: **Elastic net bandage. Roll with 25m, stretched. Size: 6 – Non Sterile**

LOT: **743738/1/A**

ARTICLE: **421-006**

SPECIFICATION	STANDARD	RESULTS
Composition	70% Polyamide 30% Polyurethane	<i>CONFORM</i>
Colour	White	<i>CONFORM</i>
Fluorescence	Absence	<i>CONFORM</i>
Length	25m ± 10%	25,2m
Width	5,5cm ± 10%	5,5cm
Weigth/ pc	>187,5g	187,8g

Approval:

Handwritten signature
Quality Management system
Sistema de gestão da qualidade
22-08-2013