

CERTIFICATE OF ANALYSIS



Issued on the	07/02/2017			
Customer	Duc Minh Long Co.			
ANALISYS RESULT FOR	EUROSILK 5 mt x 2,5 cm			
BATCH N°	17/00207			
REF.	707002			
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT	ANALISYS METHOD
Weight support	Pass	70 g/m ²	80 g/m ²	/
Colour	Pass	white	white	/
Adhesive	Pass	55 g/m ²	65 g/m ²	/
Adhesiveness180°	Pass	3.00 N/2,5 cm	5.00 N/2,5 cm	F.U IX ediz.
We certify that the above mentioned goods have been tested and they comply to the technical specifications				

EUROFARM S.p.A.
 QA DR GIOVANNI TROINA
EUROFARM. S.p.A. a Socio Unico
 Zona Industriale - Piano Tavola
 95032 BELPASSO (CT)
 P. IVA 0075 3720 879

CERTIFICATE of ANALYSIS



Issued on the	07/02/2017		
Customer	Duc Minh Long		
ANALYSIS RESULT FOR	Silvermed s.b.	cm 10 x 10 cm	
BATCH N°	16/01942		
REF.	892073		
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT
Weight support	Pass	265 g/m ²	325 g/m ²
Absorption	Pass	800%	/
Time of absorption	Pass	/	1 sec
We certify that the above mentioned goods have been tested and they comply to the technical specifications			

EUROFARM S.p.A.
 QA DR GIOVANNI TRONA
 EUROFARM. S.p.A. a Socio Unico
 Zona Industriale - Piano Tavola
 95032 Belpasso (CT)
 P. IVA 0075 3720 879

CERTIFICATE of ANALYSIS



Issued on the	07/02/2017			
Customer	Duc Minh Long Co.			
ANALYSIS RESULT FOR	Euomed cm 5 x 7 cm			
BATCH N°	17/00179			
REF.	717800			
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT	ANALYSIS METHOD
Weight support	Pass	42 g/m ²	52 g/m ²	/
Colour	white	white	white	/
Adhesive	Pass	25 g/m ²	35 g/m ²	/
Adhesiveness 180°	Pass	3.50 N/2,5 cm	5.50 N/2,5 cm	F.U.I. IX ediz.
We certify that the above mentioned goods have been tested and they comply to the technical specifications				

EUROFARM S.p.A.

QA DR GIOVANNI TROINA

EUROFARM S.p.A. a Socio Unico

Zona Industriale - Piano Tavola

95032 BELPASSO (CT)

P. IVA 0075 3720 879

CERTIFICATE OF ANALYSIS



Issued on the	07/02/2017			
Customer	Duc Minh Long			
ANALYSIS RESULT FOR	Euromed cm 12 x 7 cm			
BATCH N°	17/00143			
REF.	717804			
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT	ANALYSIS METHOD
Weight support	Pass	42 g/m ²	52 g/m ²	/
Colour	Pass	white	white	/
Adhesive	Pass	25 g/m ²	35 g/m ²	/
Adhesiveness 180°	Pass	3.50 N/2,5 cm	5.50 N/2,5 cm	F.U.I. IX ediz.
We certify that the above mentioned goods have been tested and they comply to the technical specifications				

EUROFARM S.p.A.
 QA DR GIOVANNI TROINA
 EUROFARM S.p.A. a Socio Unico
 Zona Industriale - Piano Tavola
 95032 BELPASSO (CT)
 P. IVA 0075 3720 879

CERTIFICATE OF ANALYSIS



Issued on the	07/02/2017			
Customer	Duc Minh Long			
ANALYSIS RESULT FOR	Euromed cm 15 x 8 cm			
BATCH N°	17/00046			
REF.	717805			
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT	ANALYSIS METHOD
Weight support	Pass	42 g/m ²	52 g/m ²	/
Colour	Pass	white	white	/
Adhesive	Pass	25 g/m ²	35 g/m ²	/
Adhesiveness 180°	Pass	3.50 N/2,5 cm	5.50 N/2,5 cm	F.U.I. IX ediz.
We certify that the above mentioned goods have been tested and they comply to the technical specifications				

EUROFARM S.p.A.
 QA DR GIOVANNI TROINA

EUROFARM S.p.A. - Succo Unico
 Zona Industriale - Piano Tavola
95032 BELPASSO (CT)
P. IVA 0075 3720 879

CERTIFICATE OF ANALYSIS



Issued on the	07/02/2017			
Customer	Duc Minh Long Co.			
ANALYSIS RESULT FOR	Euromed cm 10 x 25 cm			
BATCH N°	16/02338			
REF.	717807			
TEST	ANALYSIS OUTCOMES	LOWER LIMIT	HIGHER LIMIT	ANALYSIS METHOD
Weight support	Pass	42 g/m ²	52 g/m ²	/
Colour	Pass	white	white	/
Adhesive	Pass	25 g/m ²	35 g/m ²	/
Adhesiveness 180°	Pass	3.50 N/2,5 cm	5.50 N/2,5 cm	F.U.I. IX ediz.
We certify that the above mentioned goods have been tested and they comply to the technical specifications				

EUROFARM S.p.A.
QA DR GIOVANNI TROINA

EUROFARM. S.p.A. a Socio Unico
Zona Industriale · Piano Tavola
95032 Belpasso (CT)
P. IVA 0075 3720 879