




# GIMA

## STEAM TAPE AUTOCLAVE INDICATOR - 19mm x 50m

**Code:** 35846  
**Category:** Sterilization Tests  
**Unit of sale:** 12 rolls  
**Minimum order:** 1



**EAN13:** 8001814292072

**Description:**

Gima range of tapes is specially formulated for packing purposes as well as to show visibly that the pack has undergone a sterilization cycle. The backing of the different tapes is made of strong, semicreped, impregnated paper which gives excellent conformability during the application. Sterile tapes unwind easily and can be torn by hand without difficulty. Special adhesives formulated to withstand high temperature and humidity ensure excellent bond to paper, plastic, nonwovens, board, metal, glass and linen. Each tape is individually packed with due date. Made in Italy.

- For steam sterilization at 134° brown strips appear in 3 min.
- For sterilization at 121° brown strips appear in 15 min.
- MEET EN ISO 11140-1 : 2005 PART RELATED TO CLASS 1



# GIMA

## DRY HEAT TAPE 19 mm - hot air sterilizer

**Code:** 35842  
**Category:** Sterilization Tests  
**Unit of sale:** 12 rolls  
**Minimum order:** 1



**EAN13:** 8001814308575

**Description:**

Gima range of tapes is specially formulated for packing purposes as well as to show visibly that the pack has undergone a sterilization cycle. The backing of the different tapes is made of strong, semicreped, impregnated paper which gives excellent conformability during the application. Sterile tapes unwind easily and can be torn by hand without difficulty. Special adhesives formulated to withstand high temperature and humidity ensure excellent bond to paper, plastic, nonwovens, board, metal, glass and linen. Each tape is individually packed with due date. Made in Italy.

- For steam sterilization at 134° brown strips appear in 3 min.
- For sterilization at 121° brown strips appear in 15 min.
- MEET EN ISO 11140-1 : 2005 PART RELATED TO CLASS 1