



GIMA

GIMA BRUSH - conical - sterile - 500 pcs.

Code: 29742
Category: Brushes and collectors
Unit of sale: box of 500 pcs.
Minimum order: 1
Type: Medical device
Class: I S
NSIS: 40860



EAN13: 8023279297423

Description: These single-use endocervical brushes (20 cm long, 6 cm Ø brush) are an efficient method of collecting vulval, endocervical and nasal surface samples for cytological/microbiological examinations.
Made in Italy. Individually packed in sterile plastic bag; minimum order: 500 pcs.