

Tên hàng: Dây oxy đường mũi (ống thông mũi)

Chủng loại: Nasal Cannula; Nasal Oxygen Cannula; Oxygen Nasal Cannula

Số công bố TTBYT:

Chủ sở hữu số công bố, nhà nhập khẩu : Công Ty TNHH TM - DV KT Hoàng Lộc

Địa chỉ: 48 Long Hưng, Phường 7, Quận Tân Bình, TP.HCM

Điện thoại: (028) 39716873 -74 - Fax: (028) 39716873-74

Nhà SX:

1. Ningbo Greatcare Trading Co., Ltd

Unit 93, Building 12, No. 818, Qiming Road, Yinzhou 315105 Ningbo, Zhejiang China

2. Ningbo Greatcare Medical Instruments Co., Ltd.

No. 92 Cishanhe Road, Chunxiao Town, Beilun District, Ningbo, China

3. Well Lead Medical Co., Ltd

C-4 Jinhua Industrial Estate, Hualong 511434 Panyu, Guangzhou, China

4. Ningbo Shengyurui Medical Appliances Co., Ltd

No. 138, Binhaisi Road, Hangzhou Bay New Zone, Ningbo City, Zhejiang Province, 315336, P.R. China

Chủ sở hữu sản phẩm: Ningbo Greatcare Trading Co., Ltd

Unit 93, Building 12, No. 818, Qiming Road, Yinzhou 315105 Ningbo, Zhejiang China

Xuất xứ: China

Thông tin ngày SX, HSD, số lot, số seri: xem trên bao bì

Thông tin bảo hành, HSD, TLKT sửa chữa sản phẩm: vui lòng liên hệ NNK

Bảo quản: nơi khô ráo, thoáng mát, tránh ánh nắng mặt trời chiếu trực tiếp

Đại diện hợp pháp công ty



Cao Xuân Quỳnh Anh

TOP CARE

NASAL OXYGEN CANNULA


Directions For Use:

1. Attach the oxygen supply tubing to the oxygen source.
2. Set the oxygen flow as per the prescribed whit.
3. Insert the nasal tips into the nostrils passing the two plastic tubes over the ears and under the chin.

Note: To fit smaller patients use not intended for reprocessing.

CAUTION:

- * For single use. Discard after use.
- * Do not store in direct sunlight, at extreme temperatures or in high humidity. Store in dry conditions.
- * Do not use if package is open or damaged.

 Ningbo Greatcare Trading Co., Ltd.
Unit 93, Building 12, No. 818, Qiming Road,
Yinzhou 315105 Ningbo, Zhejiang, China

 Greatcare Medical GmbH
Bonner Str. 31, 50389 Wesseling, Germany

CE 0197

[LOT] XXXXX

 YYYY-MM

 YYYY-MM

SINGLE USE ONLY

[STERILE] [EO]

SIZE: L-S-XS

MADE IN CHINA