

PRODUCT NAME	Washing Solution
CATALOG NUMBER	1030165
CONTAINER SIZE	200ml of Nonionic detergent (Triton™ X-100, 0.1%)


INTENDED USE

Reagent for washing flow path of Biolyte V, Electrolyte Analyzer System.

REAGENT PREPARATION AND STABILITY

Reagent are ready for use and stable until the expiration date stated on the label, when stored at 15~30°C.

INFORMATION

Please refer to the Biolyte V REFERENCE MANUAL for application of the above reagent.



P/N : MA9010178-02 A1

