

## Digital thermometer

DT007

### FUNCTION DISPLAY

- Oral, underarm or rectal
- Fever alarm
- LCD screen
- °C / °F switch
- Beep function
- Auto off
- Cell battery

### TECHNICAL SPECIFICATIONS

1. Normal using condition:

Ambient temperature: 41°F~104°F(5°C~40°C)

Relative humidity: ≤85%

Pressure altitude: 700hPa -1060hPa

2. Storage and shipping condition:

Ambient temperature: -4°F~131°F(-20°C~55°C)

Relative humidity: ≤95%

Pressure altitude: 700hPa 1060hPa

3. Power: DC1. 5V (1pc AG3 button cell)

4. Measuring time: in 30secs(ambient temperature 77°F(25°C),  
standard water bath)

5. Measuring range: 89.6°F~109.2°F(32.0°C~42.9°C)

6. Measuring precision:

89.6°F~ 96.6°F ( 32.0°C~35.9°C ) ±0.4°F ( ±0.2°C )

96.8°F~102.2°F ( 36.0°C~39.0°C ) ±0.2°F ( ±0.1°C )

102.4°F~109.2°F ( 39.1°C~42.9°C ) ±0.4°F ( ± 0.2°C )

7. Automatic power-off : ≤10mins



**Berrcom**<sup>®</sup>

**GUANGZHOU BERRCOM MEDICAL DEVICE CO., LTD.**

Address: No.38 Huanzhen Xi Road, Dagang Town, Nansha,  
511470 Guangzhou, Guangdong, PEOPLE'S REPUBLIC OF CHINA

Tel: +86(20)34938449

Fax: +86(20)34936960



Please read the product manual carefully or purchase and use under the guidance of medical staff  
For details on contraindications or precautions