

CREJ2

cobas c 111

H314, P280, P301+P330+P331, P303+P361+P353, P304+P340, P305+P351+P338, P310, P501

REF 05401755190

CONTENT

R1 REAGENT 4 x 20 mL

SR REAGENT 4 x 5.3 mL

Creatinine Jaffé Gen.2 (compensated)
COBAS and COBAS C are trademarks of Roche.

Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN

Made in Germany

07020278001(1)

04

IVD CE 0123 15-25 °C

dialog.roche.com

Roche

cobas®

CREJ2

REF (240)05401755190

GTIN (01)04015630926824

LOT (10)00000000

UDI

2000-01-31

2000-01-01

05 Roche

For USA:

CONTENT

R1 Potassium hydroxide 900 mmol/L; phosphate 135 mmol/L, pH ≥13.5

SR Picric acid 38 mmol/L, pH 6.5

Rx only

05401755 190

002 Roche

CREJ2 1

H314

05401755

20 mL

REAGENT

15-25 °C

IVD

LOT 00000000

2000-01-31

12311234567812120000001

07021160(1)

002 Roche

2 CREJ2 2

05401755

5.3 mL

REAGENT

15-25 °C

IVD

LOT 00000000

2000-01-31

12311234567812120000001

001 Roche

**Nhãn sản phẩm
CREJ2 (05401755190)**

**Thuốc thử xét nghiệm định lượng creatinine
CREJ2**

Mã sản phẩm: 05401755190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991