



# Certificate of Analysis

Abbott Ireland Diagnostics Division,  
Lisnamuck,  
Longford,  
Co. Longford,  
Ireland

**Product Name:** Alinity i Urine NGAL Calibrators

List number: 09P2901

Lot number: 33703UD00

Expiry Date: 2023-03-15

SCP / revision

9P290001/002

Component name	List Number	Lot Number	Production date	Expiry date
Alinity i Urine NGAL CAL A	9P29A	33695UD00	2021-10-15	2023-04-15
Alinity i Urine NGAL CAL B	9P29B	33696UD00	2021-10-12	2023-04-12
Alinity i Urine NGAL CAL C	9P29C	33697UD00	2021-10-12	2023-04-12
Alinity i Urine NGAL CAL D	9P29D	33698UD00	2021-10-12	2023-04-12
Alinity i Urine NGAL CAL E	9P29E	33699UD00	2021-10-12	2023-04-12
Alinity i Urine NGAL CAL F	9P29F	33700UD00	2021-09-15	2023-03-15

## Acceptance Table: Reported Result (ng/mL)

Description	Lower Limit	Upper Limit	Value	Result
Low Control	17.00	22.89	19.52	Acceptable
Medium Control	170.00	230.00	200.10	Acceptable
High Control	1020.00	1293.00	1202.14	Acceptable

This product has been manufactured, including packaging and quality control in full compliance with ISO 13485, IVD Directive 98/79/EC and GMP guidelines as applicable.

Note: The specifications in this certificate reflect the internal Abbott release specifications. When performing the assay in a laboratory the specifications listed in the package insert apply.

Quality:

*CARMA FARZHI (Farzhi)*  
Printed name and signature

Date:

22 Nov 21