

## ISE Buffer Solution

---

**NAME:** ISE Buffer Solution

**CODE:** 204-124

**PACKAGE SIZE:** 2000mL

**APPLICATION:** Used in biochemistry analyzer, for the maintenance of instrument.

**INGREDIENT:** Sodium Hypochlorite, Sodium Hydroxide, Stabilizer

**PHYSICAL STATE:** Light Yellow Clear liquid.

**STORAGE CONDITION:** Store at 2°C~ 30°C in clean place, avoid being pressed.

**PERIOD OF VALIDITY:** 24 months

**INSTRUMENT MODEL:** Roche series biochemistry analyzer

**SAMPLE REQUIREMENT:** Blood serum

**USE METHOD:** According to the operation manual of instrument.

**WARNINGS:** 1. For in vitro diagnostic use only.

2. In case of eye or skin contact, flush affected area with copious amounts of water.
3. Induce vomiting and consult with physician in case of ingestion.
4. Dispose of waste water in compliance with local regulations.

**PRODUCTION DATE & LOT NO.:** See labels or outer package.

  
GIÁM ĐỐC  
*Lãng Đỗ Dũng*