

Mẫu nhãn

- Nhãn gốc:

Greetmed Oral preformed
Tracheal Tubes

DIRECTIONS:

- Prior to intubation, deflate the cuff completely.
- After intubation, inflate the cuff using the minimum volume of air required to provide an effective seal.
- Immediately after cuff inflation, auscultate both lung fields. If breath sounds diminished over one lung field or absent over one or both fields, adjust the tube as required.
- Endotracheal tube placement should be confirmed by viewing the position of the tube tip with a chest radiograph.

WARNING:

- Single use only.
- Sterile if package is unopened or undamaged.
- Do not resterilize.
- Do not expose to temperatures above 49°C.

CODE #: GT014-500
TYPE: WITH CUFF

SIZE:
I.D.: mm
O.D.: mm

LOT:

6 495 862 715 074 09

STERILE EO

◆ Greeting the day Greeting the future

- Nhãn phụ:

Tên hàng: Ống đặt nội khí quản

Mã SP: Xem trên bao bì

Số lưu hành:

Số LOT: Xem trên nhãn gốc [LOT], NSX: [NSX], HSD: [HSD] (năm, tháng, ngày)

Nhà sản xuất/Xuất xứ: Ningbo Greetmed Medical Instruments Co., Ltd

16F-1, Building 1, No.98 Chuangyuan Road, Hi-Tech Zone, 315042 Ningbo, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA

Chủ sở hữu số lưu hành: Công ty TNHH sản xuất – thương mại – trang thiết bị y tế Phan Anh T3-T4 Hồng Lĩnh, cư xá Bắc Hải, Phường 15, Quận 10, TP. Hồ Chí Minh, Việt Nam

Thông tin về công dụng, cảnh báo, hướng dẫn sử dụng, hướng dẫn bảo quản: Tra cứu theo số lưu hành tại <https://dmecc.moh.gov.vn/van-ban-cong-bo>