

UREAL B C

UREAL

For USA: **CONTENT**
B: TRIS 220 mmol/L, pH 8.6;
 2-oxoglutarate 73 mmol/L;
 NADH 2.5 mmol/L;
 ADP 6.5 mmol/L; urease
 (jack bean) ≥ 300 µkat/L;
 GLDH (bovine liver)
 ≥ 80 µkat/L
C: NaCl 9 %

Urea/BUN
cobas c systems
REF 05171873214
 System-ID 05 6303 9

CONTENT 1900
 2-8 °C **IVD**
 dialog.roche.com

CE 0123

Roche Diagnostics GmbH
 Sandhofer Str. 116
 D-68305 Mannheim

Distribution in USA by:
 Roche Diagnostics, Indianapolis, IN
 Made in China

LOT 00000000
 2000-01-31

COBAS and COBAS C
 are trademarks of Roche.

Roche

cobas[®] 02

04432002001(3)


Pos.1 - 05171873401 - V2




Pos.2 - 05171873402 - V1

UREAL **REF** (240)05171873214 For USA: Rx only

LOT (10)00000000
 2000-01-31
 2000-01-01

UDI 

GTIN (01)07613336161937

For USA: Store upright 

EU: +800 5505 6606
 USA: +1 800 4282336

Roche

02

Pos.3 - 05171873403 - V2

Nhãn sản phẩm
UREAL (05171873214)

**Thuốc thử xét nghiệm định lượng urea/urea
nitrogen
UREAL**

Mã sản phẩm: 05171873214

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics (Suzhou) Ltd., No. 259
Zhongyuan Road, 215028 Suzhou, Jiangsu,
Trung Quốc

Xuất xứ: Trung Quốc

Hướng dẫn bảo quản: xem  , **Số lô:** xem LOT ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991