

MẪU NHÃN


Người đại diện hợp pháp của cơ sở

Nhãn phụ:

**Thuốc thử xét nghiệm định lượng Lactate
Atellica CH Lac_2**

Số lưu hành: xxxx

Số lô: Xem  trên nhãn gốc

Hạn sử dụng: Xem  trên nhãn gốc

Bảo quản: Xem  trên nhãn gốc

Chủ sở hữu TTBYT: Xem  trên nhãn gốc.

Nước chủ sở hữu: Mỹ

Xuất xứ: Vương Quốc Anh

Chủ sở hữu số lưu hành: Công ty TNHH Thiết bị y tế Phương Đông

Địa chỉ: Tầng 04, Tòa nhà D1, ô đất CT2, Khu đô thị mới Kim Văn – Kim Lũ, Phường Đại Kim, Quận Hoàng Mai, Thành phố Hà Nội, Việt Nam.

Cảnh báo, Hướng dẫn sử dụng, Hướng dẫn bảo quản xin xem trên Hướng dẫn sử dụng tại:
www.xxxx.xxx

CTN_C3 RDX_ASSAY_Lac_2_ACH
Part Number: 11417480_01
Carton: C3 - 70M2064102
Template Number: 11417483_01

COLORS:
Black
PMS 321

DIE LINE LAYER DOES NOT PRINT

Lac_2

REF 11532568

2°C - 8°C

SIEMENS
Healthineers

Atellica® CH

Lac_2

Lac_2 R1 4 x 1.5 g

Σ 380 ㊦

IVD ⓘ siemens.com/eifu

For US / CA

RxOnly

The Atellica® CH Lactate_2 (Lac_2) Assay is for *in vitro* diagnostic use in the quantitative determination of lactate in human plasma (lithium heparin, potassium EDTA and sodium fluoride) and cerebrospinal fluid (CSF) using the Atellica® CH Analyzer.
Contains: P1 [PIPES buffer (pH 7.2) (100 mmol/L); TOOS (2.1 mmol/L); sodium azide (< 0.1%)] R1 [Peroxidase (≥ 1000 U/L); lactate oxidase (≥ 600 U/L); ascorbate oxidase (≥ 10,000 U/L); 4-aminoantipyrine (0.4 mmol/L)]

SIEMENS
Healthineers

Atellica® CH

Lac_2

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591 USA
siemens-healthineers.com

EC REP
Siemens Healthcare Diagnostics
Manufacturing Ltd.
Chapel Lane
Swords, Co. Dublin, Ireland




08E ㊦

Lac_2

Atellica® CH

SIEMENS
Healthineers

11417480_01

Placeholder 

Superbrand PMS 321
Preprinted

Product, Assay Mnemonic
Arial Narrow Bold 14 pt

Lot / Expiration Symbol
Lot Number and Expiration Date
Arial Narrow Bold, 7 pt

IUO Symbol, For Pre-Commercial Use Only

Preprinted Symbols
Single, Double, Triple

Legal Manufacturer
Preprinted


Artwork PN, Rev. Numerical (2 Digit)
Arial Narrow 6 pt


Atellica® CH

Lac_2 P1

LOT 製造番号 IVD

XXXXXX

⌚ 使用期限  95

YYYY-MM-DD  2°C 8°C

XXXXX

Siemens Healthcare Diagnostics Inc.
Tarrytown, NY 10591 USA

11417137 Rev. 01

xxxx xxxx xxxx xxxx xxxx xxx

IVD Symbol, For Commercial Use Only

Contains Sufficient XXXX Test Symbol
Numbers, Arial Narrow Bold, 7 pt

Temperature Limitation Symbol
Numbers, Arial Narrow Bold, 6 pt


Label Sequence No.
Zebra Triumvirate 6 pt

Signal Word (above or below GHS Symbols)
Arial Narrow Bold, 7 pt
WARNING or DANGER

Barcode Human Readable
Zebra Triumvirate 5 pt

AW,AT,LABEL,P1,TEST,Lac_2
 Part Number: 11417137 Rev. 01
 Label Stock Number: 11109590 Rev. 03
 Template Number: 11110965 Rev. 05

COLORS:
 Black
 PMS 321



*Lot and Expiration Date
Arial Narrow, 7 pt*

Label Specifications

Type of artwork: Label

Description: LABEL,VIAL,RX,2.67,Lac_2 R1

Part Number: 11417139 Rev. 01

Dimensions: 0.866" x 2.677"w CR .125"

Colors:

Black

PMS 321

NOTE: Key line does not print