

CLI

H314, P264, P280, P303+P361+P531
P305+P351+P338, P310, P501
EUH031

sysmex

Components:

Sodium hypochlorite 1.0%
(available chlorine concentration)

Composants:

Hypochlorite de sodium 1.0%
(Concentration du chlore a)

Componentes:

Hipoclorito de sodio 1.0%
(concentración en cloro del)

Composição:

Hipoclorito de sódio 1.0%
(concentração de cloro disponível)

組成:

次亜塩素酸ナトリウム 1.0%
(有効塩素濃度)

GSA-500A
54212541G

(体外診断用試剤)
全自動血液凝固分析装置用洗液
危険

造成严重皮肤灼伤和眼损伤
应急咨询电话: 400-820-0815
警告及注意事项:
使用前请仔细阅读该包装说明书及化学品安全技术说明书。

主要组成成分:
次亜塩素酸 (指鉄濃度) 1.0%

经销商及售后服务单位:
希森美康医疗电子(上海)有限公司
上海市浦东新区川沙路1233号汇亚大厦9楼(200120)
电话: 400-820-0815
生产商: SYSMEX CORPORATION
注册地址: 1-5-1, Wakinohama-Kaigandori,
Chuo-ku, Kobe, 651-0073, Japan
电话: 81-78-265-0500
生产地址: 17 Takumidai, Ono, Hyogo, 675-1322,
Japan

sysmex

CLI CA CLEAN® I
CAクリーン® I

REF 964-0631-3 50mLx1

IVD +2°C ~ +8°C CE

i www.sysmex-europe.com/ifu

EC REP Sysmex Europe GmbH
Bainbarch 1,
22848 Norderstedt, Germany

Sysmex Corporation
1-5-1 Wakinohama-Kaigandori,
Chuo-ku, Kobe 651-0073, Japan
MADE IN JAPAN

Detergent for fully automated blood coagulation analyzer
Détergent pour analyseur de coagulation sanguine entièrement automatisé
Detergente para analizador de coagulación de sangre totalmente automático
Detergente para analisador totalmente automático de coagulação do sangue
全自動血液凝固測定装置用洗液

Distributed in the U.S.A. by
Sysmex America, Inc.
577 Aptakissic Road, Lincolnshire, IL 60069, U.S.A.
Tel: +1-224-543-9500
Sysmex is Registered in U.S. Pat. & Tm. Off.

GHS requirements:
Danger
H314 Causes severe skin burns and eye damage.
P264 Wash skin thoroughly after handling.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P303+P361+P531 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.
P501 Dispose of contents/container in accordance with local and national regulations.

危険
重篤な皮膚の灼傷・眼の損傷。
取扱後は、皮膚をよく洗うこと。
保護手袋/保護衣/保護服/保護面を着用すること。
皮膚(または髪)に付着した場合、着ちに汚染された衣服をすべて脱ぐこと。皮膚を流水/シャワーで洗うこと。
顔に入った場合、水で数分間注目を流すこと。次に、コンタクトレンズを適用して、容易に外せる場合は除外すること。その他も洗浄を続けること。
直ちに医師に連絡すること。
内容物/容器を都道府県知事の許可を受けた専門の廃棄処理業者に委託して、廃棄すること。
Before use, please read the package insert and safety data sheet carefully.
Avant utilisation, se reporter aux informations de la notice et aux fiches de données de sécurité.
Antes de usar, leer atentamente las instrucciones del interior y ficha de datos de seguridad.
Antes de usar, ler atentamente a bula e ficha de segurança.
ご使用の前に、添付文書及び安全データシートを必ずお読みください。
製造販売元
シスメックス株式会社
神戸市中央区臨海町岸部1丁目5番1号
Tel: 078-265-0500



GIÁM ĐỐC
Lê Trung Hữu

印刷色	
	PANTONE 376C
	PANTONE 1795C
	DIC 139
	黒

CA CLEAN I_GSA-500A_54212541G_404306033_Box