






Labels






A. Outer Carton Box






Product	Medical Head Lamp	Model	L2SW
Rating	3.7 V  , 350 mA	SN	mark on the inner box
Weight	kg		
Quantity	1 set	EC REP	Cormedics Medizintechnik GmbH Bahnhofstr. 32 82041 Delsenhofen, Germany TEL. +49 (0)89 666 58 79 15 Fax. +49 (0)89 54 999 030
	XENOSYS CO., LTD 2, Venture-ro 100beon-gil, Yeonsu-gu, Incheon 22013, Korea TEL. +82-32-875-9811		
The medical device			  CE

B. LED Control Unit







Model LCUW (Control Unit)
 Rating 3.7 V , 350 mA
 Xenosys Co., Ltd.
 www.xenosys.co.kr   CE
 Tel. +82-32-875-9811
 SN

C. Battery Pack

 Model BTPWL (1ICP7/20/45)
 Rating 3.7 V , 620 mAh (2.29 Wh)
 www.xenosys.co.kr   CE
 SN Secondary LI

 Model BTPWS (1ICP7/20/30)
 Rating 3.7 V , 380 mAh (1.41 Wh)
 www.xenosys.co.kr   CE
 SN Secondary LI

D. Charger

Product	Medical Head Lamp(LCHW)	Model	L2SW
Certl. No.	Gyeongin jesin 07-86	Quantity	1 SET
SN			
Power	Input: 5 V  , 2 A Output: 4.2 V  , 0.70 A	The medical device	
	XENOSYS CO., LTD 2, Venture-ro 100beon-gil, Yeonsu-gu, Incheon 22013, Korea TEL. +82-32-875-9811	  CE	