



Labels

A. Outer Carton Box

Product	Medical Head Lamp	Model	L2SWE
Rating	3.7 V , 350 mA	SN	mark on the inner box
Weight	kg		
Quantity	1 set	EC REP	Cormedics Medzintechnik GmbH Bahnhofstr. 32 82041 Delsenhofen, Germany TEL. +49 (0)89 666 58 79 15 Fax. +49 (0)89 54 999 030
	XENOSYS CO., LTD 2, Venture-ro 100beon-gil, Yeonsu-gu, Incheon 22013, Korea TEL. +82-32-875-9811		
The medical device			CE

B. LED Control Unit

Model LCUW (Control Unit)
 Rating 3.7 V , 350 mA

Xenosys Co., Ltd. CE
 www.xenosys.co.kr
 Tel. +82-32-875-9811

SN

C. Battery Pack

Model BTPWS (652030 LI polymer)
 Rating 3.7 V , 380 mAh (1.41 Wh)
 www.xenosys.co.kr CE

SN

Model BTPWL (652045 LI polymer)
 Rating 3.7 V , 620 mAh (2.29 Wh)
 www.xenosys.co.kr CE

SN

D. Charger

Product	Medical Head Lamp(LCHW)	Model	L2SWE
Certl. No.	Gyeongin jesin 07-86	Quantity	1 SET
SN			
Power	Input: 5 V , 2 A Output: 4.2 V , 0.70 A	The medical device	
	XENOSYS CO., LTD 2, Venture-ro 100beon-gil, Yeonsu-gu, Incheon 22013, Korea TEL. +82-32-875-9811		CE