

Yongyue Medical Technology(Kunshan) Co.,Ltd.

Certificate of Compliance

Product Name:10ul pipette tip

Cat.No.:	YPT0010-T-F-SL	Lot No.:	E210729-0102
Certified:	ISO13485:2016		

Description

Material:	Polypropylene(Metal free)	Pass
Request:	DNA/Proteins free	Pass

Sampling item	Criterion	Result
	Sampling rate:1%.	
Appearance	Requirements:No defect,no bubbles,no burrs,no pollution,no deformation.	Pass
Dimensions	Length:38.95-39.05 mm,width:3.23-3.33mm	Pass

We guarantee that this product is produced in a class 6(according to VDI 2083)and a class 100,000(according to U.S.Fed.Stand.209D) clean room environment,and fulfils the following purity.

Specification	Method	Test result
RNase	RNA digestion	< 1.0x10 ⁻⁶ Kunitz units
DNase	DNA digestion	< 1.0x10 ⁻⁶ Kunitz units
Pyrogen	Polymerase Chain Reaction	<=0.1 EU/ml

Final result	Pass
---------------------	-------------

QC: Liusihao
Data: 2021.07.29

Approved: *Witt.xie*
Data:2021.07.29



Tel:+86-512-55256601-8870

Fax:+86-512-57260163

E-mail:qs@kszyg.com

No. 6 Kingdee Road, Yushan Town, Kunshan, Suzhou City, Jiangsu province, P. R. China



Yongyue Medical Technology(Kunshan) Co.,Ltd.
Certificate of Compliance

Product Name:20ul pipette tip

Cat.No.: YPT0020-T-F-SL Lot No.: E210729-0202
Certified: ISO13485:2016

Description

Material: Polypropylene(Metal free) Pass
Request: DNA/Proteins free Pass

Sampling item	Criterion	Result
	Sampling rate:1%.	
Appearance	Requirements:No defect,no bubbles,no burrs,no pollution,no deformation.	Pass
Dimensions	Length:55.95-56.05mm,width:5.22-5.32mm	Pass

We guarantee that this product is produced in a class 6(according to VDI 2083)and a class 100,000(according to U.S.Fed.Stand.209D) clean room environment,and fulfils the following purity.

Specification	Method	Test result
RNase	RNA digestion	<1.0x10 ⁻⁶ Kunitz units
DNase	DNA digestion	<1.0x10 ⁻⁶ Kunitz units
Pyrogen	Polymerase Chain Reaction	<=0.1 EU/ml

Final result **Pass**

QC: Liusihao
Data: 2021.07.29

Approved: *Witt.xie*
Data:2021.07.29



Tel:+86-512-55256601-8870

Fax:+86-512-57260163

E-mail: @kszy.com

No. 6 Kingdee Road, Yushan Town, Kunshan, Suzhou City,
Jiangsu province, P. R. China



Yongyue Medical Technology(Kunshan) Co.,Ltd.

Certificate of Compliance

Product Name:100ul pipette tip

Cat.No.:	YPT0100-T-F-SL	Lot No.:	E210729-0302
Certified:	ISO13485:2016		

Description

Material:	Polypropylene(Metal free)	Pass
Request:	DNA/Proteins free	Pass

Sampling item	Criterion	Result
	Sampling rate:1%.	
Appearance	Requirements:No defect,no bubbles,no burrs,no pollution,no deformation.	Pass
Dimensions	Length:55.95-56.05mm,width:5.22-5.32mm	Pass

We guarantee that this product is produced in a class 6(according to VDI 2083)and a class 100,000(according to U.S.Fed.Stand.209D) clean room environment,and fulfils the following purity.

Specification	Method	Test result
RNase	RNA digestion	< 1.0x10 ⁻⁶ Kunitz units
DNase	DNA digestion	< 1.0x10 ⁻⁶ Kunitz units
Pyrogen	Polymerase Chain Reaction	<=0.1 EU/ml

Final result	Pass
---------------------	-------------

QC: Liusihao
Data: 2021.07.29

Approved: *Witt. xie*
Data:2021.07.29



Tel:+86-512-55256601-8870

Fax:+86-512-57260163

E-mail:qs@kszyg.com

No. 6 Kingdee Road, Yushan Town, Kunshan, Suzhou City, Jiangsu province, P. R. China



Yongyue Medical Technology(Kunshan) Co.,Ltd.

Certificate of Compliance

Product Name:200ul pipette tip

Cat.No.:	YPT0200-T-F-SL	Lot No.:	E210729-0402
Certified:	ISO13485:2016		

Description

Material:	Polypropylene(Metal free)	Pass
Request:	DNA/Proteins free	Pass

Sampling item	Criterion	Result
	Sampling rate:1%.	
Appearance	Requirements:No defect,no bubbles,no burrs,no pollution,no deformation.	Pass
Dimensions	Length:55.95-56.05mm,width:5.22-5.32mm	Pass

We guarantee that this product is produced in a class 6(according to VDI 2083)and a class 100,000(according to U.S.Fed.Stand.209D) clean room environment,and fulfils the following purity.

Specification	Method	Test result
RNase	RNA digestion	< 1.0x10 ⁻⁶ Kunitz units
DNase	DNA digestion	< 1.0x10 ⁻⁶ Kunitz units
Pyrogen	Polymerase Chain Reaction	<=0.1 EU/ml

Final result	Pass
---------------------	-------------

QC: Liusihao
Data: 2021.07.29

Approved: *Witt.xie*
Data:2021.07.29



Tel:+86-512-55256601-8870

Fax +86-512-57260163

E-mail: cs@kszy.com

No. 6 Kingdee Road, Yushan Town, Kunshan, Suzhou City, Jiangsu province, P. R. China



Yongyue Medical Technology(Kunshan) Co.,Ltd.

Certificate of Compliance

Product Name:1000ul pipette tip

Cat.No.:	YPT1000-T-F-SL	Lot No.:	E210729-0502
Certified:	ISO13485:2016		

Description

Material:	Polypropylene(Metal free)	Pass
Request:	DNA/Proteins free	Pass

Sampling item	Criterion	Result
	Sampling rate:1%.	
Appearance	Requirements:No defect,no bubbles,no burrs,no pollution,no deformation.	Pass
Dimensions	Length: 99.95-100.05mm,width: 7.85-7.89mm	Pass

We guarantee that this product is produced in a class 6(according to VDI 2083)and a class 100,000(according to U.S.Fed.Stand.209D) clean room environment,and fulfils the following purity.

Specification	Method	Test result
RNase	RNA digestion	< 1.0x10 ⁻⁶ Kunitz units
DNase	DNA digestion	< 1.0x10 ⁻⁶ Kunitz units
Pyrogen	Polymerase Chain Reaction	<=0.1 EU/ml

Final result	Pass
---------------------	-------------

QC: Liusihao
Data: 2021.07.29

Approved: *Witt.xie*
Data:2021.07.29



Tel:+86-512-55256601-8870

Fax:+86-512-57260163

E-mail: cs@yongyue.com

No. 6 Kingdee Road, Yushan Town, Kunshan, Suzhou City, Jiangsu province, P. R. China

