

ALB2

B
ALB2

C

For USA: **CONTENT**

B: Citrate buffer
95 mmol/L, pH 4.1

C: Citrate buffer
95 mmol/L, pH 4.1;
bromocresol green
0.66 mmol/L

Albumin Gen.2
cobas c systems

REF 08056692190

System-ID 20 0900 1

CONTENT 750

15-25 °C **IVD**

 dialog.roche.com

COBAS and COBAS C
are trademarks of Roche.



 Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
Made in Germany

cobas[®]

04432002001(3)

LOT 00000000

 2000-01-31



ALB2

REF (240)08056692190

For USA: Rx only

LOT (10)00000000

 2000-01-31


 2000-01-01



For USA: Store upright

UDI



 EU: +800 5505 6606
USA: +1 800 4282336

GTIN (01)07613336120941



03

**Nhãn sản phẩm
ALB2 (08056692190)**

**Thuốc thử xét nghiệm định lượng albumin
ALB2**

Mã sản phẩm: 08056692190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991