

ALB2

B
ALB2

C

For USA: **CONTENT**
B: Citrate buffer 95 mmol/L,
pH 4.1
C: Citrate buffer 95 mmol/L,
pH 4.1; bromocresol green
0.66 mmol/L

Albumin Gen.2
COBAS INTEGRA / **cobas c** systems

REF 03183688122
System-ID 07 6592 9

CONTENT 300

15-25 °C **IVD**

 dialog.roche.com

COBAS INTEGRA, COBAS C and
COBAS are trademarks of Roche.



Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
Made in Germany


04432002001(3)

LOT 00000000

 2000-01-31



cobas[®]




Ser.No. 00001 **ALB2** 03183688 122 07

ALB2 **REF** (240)03183688122 For USA: Rx only

 2000-01-31

 2000-01-01

 For USA: Store upright

UDI



 EU: +800 5505 6606
USA: +1 800 4282336

LOT (10)00000000

GTIN (01)04015630918898



10

**Nhãn sản phẩm
ALB2 (03183688122)**

**Thuốc thử xét nghiệm định lượng albumin
ALB2**

Mã sản phẩm: 03183688122

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,


Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991