



## Document Information:

Document Author: Scholl, Alexandra {DOMC~Mannheim}  
Business Area / Unit: Roche Professional Diagnostics  
Confidentiality: Confidential  
Document Class: Device Master Record  
Document Type: Device Master Record  
Document Creator: Scholl, Alexandra {DOMC~Mannheim}  
Document Lifecycle Status: Signed  
Valid From: 08-Oct-2020 15:00:29 (UTC)  
Valid To:  
Document Title: C\_0201\_RSP\_08056692000  
Document Number: 0000000000001200000303639  
Document Version: 02  
Template: No

Global Group: Regulatory Affairs  
Global SubGroup: Specification Prod. (Registr.)  
Local Group: Reg. Specification  
Site: RDG Germany  
Department: C\_DOMC  
Document Description: Albumin Gen.2 (ALB2) - "Only for use in Regulatory Affairs"  
Change Request eCCP PLM: 400133569

## Electronic Signatures:

Signed By: scholla (Alexandra Scholl {DOMCSD})  
Role: Reviewer  
Signature Differentiation: Quality Control  
Signed Date: 08-Oct-2020 11:38:55 (UTC)

Signed By: kneisses (Sandra Kneissel {DOMC})  
Role: Approver  
Signature Differentiation: Quality Control  
Signed Date: 08-Oct-2020 14:55:53 (UTC)



## Specification

Cat. No.: 08056692190

Product: **Albumin Gen.2 (ALB2)**

<b>Parameter</b>	<b>Specification</b>
<b>Reagents-working solutions</b>	
<b>Buffer solution (B):</b>	
pH-value (25 °C):	4.00 - 4.20
<b>BCG/buffer solution (C):</b>	
pH-value (25 °C):	4.00 - 4.20

**Parameter****Specification**

---

**Performance test****Recovery on cobas c 303 analyzer:**

PreciControl ClinChem Multi 1 (Cat. No. e.g. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. e.g. 05947774190): 90 - 110 %

**Recovery on cobas c 503 analyzer:**

PreciControl ClinChem Multi 1 (Cat. No. e.g. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. e.g. 05947774190): 90 - 110 %

**Also required:**

Calibrator f.a.s. (Cat. No. 10759350190)

**Roche Diagnostics GmbH  
Quality Control Mannheim**Sandhofer Strasse 116  
D-68305 Mannheim  
Phone +49-621-759-0