



## Document Information:

Document Author: Maschat, Cordula {DOMC~Mannheim}  
Business Area / Unit: Roche Professional Diagnostics  
Confidentiality: Confidential  
Document Class: Device Master Record  
Document Type: Device Master Record  
Document Creator: Maschat, Cordula {DOMC~Mannheim}  
Document Lifecycle Status: Signed  
Valid From: 21-Sep-2020 11:02:46 (UTC)  
Valid To:  
Document Title: C\_0201\_RSP\_08056960000  
Document Number: 0000000000001200000306245  
Document Version: 02  
Template: No

Global Group: Regulatory Affairs  
Global SubGroup: Specification Prod. (Registr.)  
Site: RDG Germany  
Department: C\_DOMCC  
Document Description: Bilirubin Total Gen.3 (BILT3) - "only for use in Regulatory  
Affairs"  
Change Request eCCP PLM: 400133569

## Electronic Signatures:

Signed By: schnabeh (Helena Schnabel {DOMCSD})  
Role: Reviewer  
Signature Differentiation: Quality Control  
Signed Date: 21-Sep-2020 10:55:55 (UTC)

Signed By: byrnesh (Helga Byrnes {DOMC})  
Role: Approver  
Signature Differentiation: Quality Control  
Signed Date: 21-Sep-2020 10:57:05 (UTC)



## Specification

Cat. No.: 08056960190

Product: **Bilirubin Total Gen.3 (BILT3)**

### Parameter

### Specification

---

#### Reagents-working solutions

##### Detergent / buffer solution (B):

Appearance:

corresponds  
(clear, slightly yellowish  
solution)

##### 3,5-dichlorophenyl diazonium salt solution (C):

Appearance:

corresponds  
(colorless to yellow colored,  
clear solution)

<b>Parameter</b>	<b>Specification</b>
------------------	----------------------

---

**Performance test**

**Recovery on cobas c 303 analyzer:**

PreciControl ClinChem Multi 1 (Cat. No. e.g. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. e.g. 05947774190): 90 - 110 %

**Recovery on cobas c 503 analyzer:**

PreciControl ClinChem Multi 1 (Cat. No. e.g. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. e.g. 05947774190): 90 - 110 %

**Also required:**

Calibrator f.a.s. (Cat. No. 10759350190)

**Roche Diagnostics GmbH  
Quality Control Mannheim**

Sandhofer Strasse 116  
D-68305 Mannheim  
Phone +49-621-759-0