



Certificate of Analysis

Abbott Ireland Diagnostics Division,
Lisnamuck,
Longford,
Co. Longford,
Ireland

Product Name: Alinity i 25-OH Vitamin D Calibrators
List number: 08P4501
Lot number: 37548UD00
Expiry Date: 2023-04-27

SCP / revision

8P450001/003

Component name	List Number	Lot Number	Production date	Expiry date
Alinity i Vit D 25OH CAL A	8P45A	37542UD00	2022-02-09	2023-05-09
Alinity i Vit D 25OH CAL B	8P45B	37543UD00	2022-01-31	2023-04-30
Alinity i Vit D 25OH CAL C	8P45C	37544UD00	2022-01-31	2023-04-30
Alinity i Vit D 25OH CAL D	8P45D	37545UD00	2022-01-27	2023-04-27
Alinity i Vit D 25OH CAL E	8P45E	37546UD00	2022-01-27	2023-04-27
Alinity i Vit D 25OH CAL F	8P45F	37547UD00	2022-01-27	2023-04-27

Acceptance Table: Reported Result (ng/mL) (SMP Calibration Curve)

Description	Lower Limit	Upper Limit	Value	Result
Low Control	18.0	22.0	19.5	Acceptable
Medium Control	36.0	44.0	39.7	Acceptable
High Control	67.5	82.5	74.9	Acceptable

This product has been manufactured, including packaging and quality control in full compliance with ISO 13485, IVD Directive 98/79/EC and GMP guidelines as applicable.

Note: The specifications in this certificate reflect the internal Abbott release specifications. When performing the assay in a laboratory the specifications listed in the package insert apply.

Quality:

Date:

CARMA FREZTI C FEOUEN
Printed name and signature

09 March 22