



as human for human

## PARAMOUNT BED

### WHEEL CHAIR, MODEL KK-T551MB



- DIMENSIONS: 630(W) x 975(D) x Front:400 • Back:400(SH) mm
- Wheel: 22-inch
- Caster: 6-inch
- Foldable type
- Frame: Made of Aluminum
- Seat: Vinyl leather, Urethane foam (Cushion)
- Brake provided

**SAO Y BẢN CHÍNH**  
Ngày.....tháng.....năm.....



**YASUYUKI TAJIMA**  
Director