




Oxygen Concentrator

KSOC-5 Series 5L/10L/Dual flow/High pressure



Oil free Compressor



24 Hours continuous working



Wise alarm function



High purity oxygen



Three-channel oxygen filter



Optional Oxygen purity alarm

- America PSA technology offer nature oxygen
- High efficiency molecular sieve bed
- Reliable and durable oil free compressor
- Available 24 hours continuous working
- Self diagnostic system with error code indication
- Optional dual flow and flow splitter for multi people use

FUNCTIONS:

- Power off alarm, Overload protection
- High/Low pressure alarm
- Temperature alarm, Error code indication
- Optional Nebulizer, Oxygen purity alarm

SPECIFICATION:



Dual flow

Model	Oxygen purity	Flow Range	Max Noise	Input Voltage	Output Voltage	Power(a)	G.W	Size(mm)	Carton size (mm)
KSOC-5	93% ± 3%	0-5L	≤54dB	220V - 240V	40 - 70kPa	≤480W	18,3kgs	380x320x590	465x375x655