

ALP2S B C


ALP2S

For USA: **CONTENT**
B: 2-amino-2-methyl-1-propanol 1.724 mol/L, pH 10.44;
Mg acetate 3.83 mmol/L;
zinc sulfate 0.766 mmol/L;
N-(2-hydroxyethyl)-ethylene-diamine triacetic acid 3.83 mmol/L
C: p-nitrophenyl phosphate 132.8 mmol/L, pH 8.50

ALP IFCC Gen.2 Small
COBAS INTEGRA / **cobas c** systems H315,H319,P264,
P280,P302+P352,
P332+P313,P337+
P313,P362+P364

REF 03333752190
System-ID 07 6761 1

CONTENT 200
2-8°C **IVD**

 dialog.roche.com **CE** 0123

COBAS INTEGRA, COBAS C and COBAS are trademarks of Roche.


Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
Made in Germany

LOT 00000000
2000-01-31

13
06781411001(01)
cobas



Pos.1 - 03333752101 - V13





Ser.No. 00001 **ALP2S** 03333752 190 07


Pos.2 - 03333752102 - V7

ALP2S **REF** (240)03333752190 For USA: Rx only

 2000-01-31
 2000-01-01

 For USA: Store upright

UDI 

 EU: +800 5505 6606 **LOT** (10)00000000
USA: +1 800 4282336 **GTIN** (01)04015630919512

Roche 10

Pos.3 - 03333752103 - V10

**Nhãn sản phẩm
ALP2S (03333752190)**

**Thuốc thử xét nghiệm định lượng alkaline
phosphatase
ALP2S**

Mã sản phẩm: 03333752190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:

Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991