

Certificate of Analysis

Product Name: ACTIVATED ALANINE AMINOTRANSFERASE Abbott List Number: 08L9222
 OEM Part Number: 08L9222 Abbott Commodity Number: 08L9222
 OEM Lot Number: 97993UN22 Expiration Date: 2023-01-25
 Storage: 2 - 8°C Package Insert No./Rev. G95938 / R03
 Qty Manufactured: 70 Kits

Specification	Specification Limit	Result
Barcode Scan	Acceptable verification on c8000	Complies Pass
Bio-Rad L1	24 - 28 U/L	24 Pass
Bio-Rad L1	Lot Number	55911
Bio-Rad L2	94 - 110 U/L	98 Pass
Bio-Rad L2	Lot Number	55912
High Standard % Difference	-10 - 10 %	0 Pass
Micro: Beg Sample (R1)	<30 CFU/mL	0 Pass
Micro: End Sample (R1)	<30 CFU/mL	0 Pass
Micro: Beg Sample (R2)	<30 CFU/mL	0 Pass
Micro: End Sample (R2)	<30 CFU/mL	0 Pass

For in-vitro diagnostic use.

This lot meets all current release specifications per quality document MS.MASTER342.

The above lot has been manufactured in compliance with Quality Systems Regulations, passed all Thermo Fisher final release specifications, and has met all Abbott acceptance requirements.

This product contains R1 and R2 components. The component name, lot number, and expiration date matches the kit name, lot number, and expiration date listed on this COA.



Quality Representative

Date: 2022-03-18