

**Certificate of Analysis**

Product Name: ACTIVATED ALANINE AMINOTRANSFERASE	Abbott List Number: 08L9242
OEM Part Number: 08L9242	Abbott Commodity Number: 08L9242
OEM Lot Number: 97995UN22	Expiration Date: 2023-01-25
Storage: 2 - 8°C	Package Insert No./Rev. G95938 / R03
Qty Manufactured: 31 Kits	

Specification	Specification Limit	Result
Barcode Scan	Acceptable verification on c8000	Complies Pass
Bio-Rad Control L1	24 - 28 U/L	24 Pass
Bio-Rad Control L1	Lot Number	55911
Bio-Rad Control L2	94 - 110 U/L	98 Pass
Bio-Rad Control L2	Lot Number	55912
High Standard % Difference	-10 - 10 %	0 Pass
Micro: Beg Sample (R1)	<30 CFU/mL	0 Pass
Micro: End Sample (R1)	<30 CFU/mL	0 Pass
Micro: Beg Sample (R2)	<30 CFU/mL	0 Pass
Micro: End Sample (R2)	<30 CFU/mL	0 Pass

For in-vitro diagnostic use.

This lot meets all current release specifications per quality document MS.MASTER342.

The above product has been manufactured in compliance with Quality Systems Regulations, passed all Thermo Fisher final release specifications, and has met all Abbott acceptance requirements.

This product contains R1 and R2 components. The component name, lot number, and expiration date matches the kit name, lot number, and expiration date listed on this COA.



Quality Representative

Date: 2022-03-18