

Certificate of Analysis

Product Name: TRIGLYCERIDES CLINICAL CHEMISTRY	Abbott List Number: 07D7422
OEM Part Number: 07D7422	Abbott Commodity Number: 07D7422
OEM Lot Number: 38825UN22	Expiration Date: 2023-09-08
Storage: 2 - 8°C	Package Insert No./Rev. G95961 / R02
Qty Manufactured: 1,241 Kits	

Specification	Specification Limit	Result
Barcode Scan	Acceptable verification on c8000	Complies Pass
Bio-Rad L1	161 - 197 mg/dL	177 Pass
Bio-Rad L1	Lot Number	55911
Bio-Rad L2	85 - 104 mg/dL	91 Pass
Bio-Rad L2	Lot Number	55912
Standard Level 3 % Difference	-10 - 10 %	-1 Pass
Micro: Beg Sample	<30	0 Pass
Micro: End Sample	<30	1 Pass
Date of Manufacture	Report	2022-03-08

For in-vitro diagnostic use.

This product contains animal source material. The bovine material used in the manufacture of this product is derived from bovine blood collected from U.S. sourced cattle slaughtered at a USDA licensed establishment.

This lot meets all current release specifications per quality document MS.MASTER342.

The above product has been manufactured in compliance with Quality Systems Regulations, passed all Thermo Fisher final release specifications, and has met all Abbott acceptance requirements.

This product contains a single component. The component name, lot number, and expiration date matches the kit name, lot number, and expiration date listed on this COA.



Quality Representative

Date: 2022-04-01