

CÔNG TY CỔ PHẦN THIÊN TRƯỜNG

CỘNG HÒA XÃ HỘI CHỦ NGHĨA VIỆT NAM  
Độc lập - Tự do - Hạnh phúc

Số: 08/PLTBYT/TTC-IMC

Hà Nội, ngày 08 tháng 9 năm 2017

## BẢN PHÂN LOẠI TRANG THIẾT BỊ Y TẾ

Kính gửi: Công ty TNHH Thiết bị và Vật tư Y tế Đông Dương

Căn cứ Nghị định số 36/2016/NĐ-CP ngày 15 tháng 5 năm 2016 của Chính phủ về quản lý trang thiết bị y tế;

Nguyên tắc được sử dụng để phân loại: Quy tắc 6- Quy tắc phân loại trang thiết bị y tế không phải trang thiết bị chẩn đoán In vitro

Chúng tôi phân loại trang thiết bị y tế như sau:

TT	Tên trang thiết bị y tế	Chủng loại/ mã sản phẩm	Hãng/ nước sản xuất	Hãng/ nước chủ sở hữu	Phân loại
1	Dụng cụ panh, kim, kẹp dùng trong phẫu thuật mắt	Theo <i>Phụ lục phân loại trang thiết bị y tế</i>	Indo Webal Surgicals Pvt.Ltd., Ấn Độ	Indo Webal Surgicals Pvt.Ltd., Ấn Độ	A

Đại diện hợp pháp của cơ sở phân loại

Nơi nhận:

- Như kính gửi
- Lưu Cty CP Thiên Trường



GIÁM ĐỐC

*Dinh Quang Khải*

**PHỤ LỤC PHÂN LOẠI TRANG THIẾT BỊ Y TẾ**  
*SỐ 08/PLTBYT/TTC-IMC, ngày 08 tháng 9 năm 2017*

<b>TT</b>	<b>Chủng loại</b>	<b>Mã sản phẩm</b>
1	Jaffe Capsule Forcep round grasping tips extremely thin 11 mm long shanks.	IW-883
2	Polack Corneal Suturing Forceps very delicate 1mm spread 1x2 teeth 0.12mm.	IW-884
3	Harms-colibri suturing forceps very delicate 1x2 teeth, 0.12mm with tying platform.	IW-885
4	Gill-Colibri forceps very fine pointed tips with tying platform.	IW-886
5	Colibri forceps very delicate 1x2 teeth 0.12mm	IW-887
6	Botvin Iris forceps for peripheral iridectomy serrated tips 1x2 teeth	IW-888
7	Castroviejo Colibri forceps very delicate 1x2 teeth, 0.12mm	IW-889
8	Barraquer colibri-very delicate forceps 0.12mm-1x2 teeth, with tying platform 7mm long	IW-890
9	Castroviejo suturing forceps 1x2 teeth with typing platform	IW-891
10	Bishop harman tissue forceps 1x2 teeth delicate	IW-892
11	Pierse type micro forceps straight skeleton fine Available in: 18-0.5, 19-0.2, 20-0.8mm tip	IW-893
12	Bonn forceps 1x2 tip, 0.12mm with typing platform	IW-894
13	Hoskin-1 Tissue forceps fine grooved tips, straight & curved	IW-895
14	Adson tissue forceps, 1x2 teeth straight, 125mm long delicate	IW-896
15	St.Martin corneal suture forceps, straight shaft 1x2 teeth with 6.0mm long typing platform	IW-897
16	Dressing forceps serrated, delicate, straight & curved	IW-898
17	Ernest Nucleus Cracker extra delicate cross action serrated paddle shaped jaws for cracking the nucleus through a 2.5 mm incision	IW-899
18	Dodick Nucleus Cracker cross action longitudinally serrated paddle shaped jaws for cracking the nucleus through a 3.0 mm incision.	IW-900
19	Kansas Nucleus Removal Forceps extra delicate multiple teeth for extracting nucleus fragments for manual small incision surgery	IW-901

20	Livernois Iol Folding Forceps concave jaw surfaces for soft iols	IW-902
21	Steinert IOL inserting Forceps conically shaped jaws with flat inside surfaces, used for folding and inserting soft iols	IW-903
22	Lims Corneal Suture Forceps curved shaped 1x2 teeth with 5.0 Mm long tying platform Available in : 0.2mm, 0.12mm extra fine teeth	IW-905
23	Mc pherson Corneal Forceps 1x2 teeth 0.4 Mm with tying platform, big straight	IW-906
24	Mc pherson Corneal Forceps 1x2 teeth 0.4 Mm with tying platform, small straight	IW-906A
25	Mc pherson Corneal Forceps 1x2 teeth 0.4 Mm with tying platform, angled	IW-906B
26	Hess Iris Forceps 0.2 mm tip standard model for peripheral iridectomy	IW-907
27	Dastoor Superior Rectus Forceps 1x2 teeth	IW-908
28	Indo Webal Fixation Forceps straight, 1x2 & 2x3 teeth	IW-909
29	Graefe Fixation Forceps 4.5 Mm wide jaws with fine teeth	IW-910
30	Green Fixation Forceps multiple teeth , 5.0mm & 10.0 mm wide jaws	IW-911
31	Scleral Fixation / S.I.C.S. Forceps 1x2 teeth 6 mm platform curved	IW-912
32	Indo Webal Fixation Forceps 2x3 teeth, straight / angular	IW-913
33	Wills Hospital Utility Forceps with straight shaft crisscross serrated platform 13.0mm long	IW-915
34	Bonaccolto Utility Forceps with straight shaft lengthwise serrated jaws 16.0 mm long 1.7 mm wide with cross serrated tips	IW-916
35	Moorfield Conjunctival Suture Forceps straight shaft plain serrated jaws with groove	IW-917
36	Barraquer Cilia Forceps lightly angled textured jaws provide superior grasping capabilities	IW-918
37	Beer Cilia Forceps with 4.5 mm long rounded platform	IW-919
38	Arruga Capsule Forceps 2.5 mm cup curved shaped.	IW-920
39	Prince Advancement Forceps 10.0mm long jaws with three teeth and sliding lock, right / left	IW-921
40	Jameson Muscle Forceps with slide lock four teeth child size, right / left	IW-922
41	Worth Advancement Forceps 11.0mm longitudinal serrated jaws with sliding lock, right / left	IW-923

42	Berke Ptosis Clamp multy curved shafts with straight blades and longitudinal serrations with sliding lock Available in: 20.0 mm, 27.0 mm	IW-924
43	Snellen Entropiom Forceps with solid lower plate and open upper plate with locking thumb screw, right / left	IW-925
44	Ayer Chalazion Forceps with 10mm open upper plate & solid lower with locking thumb screw	IW-926
45	Lambert Chalazion Forceps round 8mm I.D	IW-927
46	Desmarres Chalazion Forceps oval 20 x13mm I.D large	IW-928
47	Hunt Chalazion Forceps round 12mm I.D	IW-929
48	Indo Webal Universal Lens Folding forceps to fold the lens perfectly and groove provided in the jaws to hold the lens in the center	IW-930
49	Shepard Lens Holding Forceps with serrated lower jaws gently curved shaft	IW-931
50	Shepard Lens Holding Forceps with 3mm lower jaw ,1.5mm upper jaw	IW-932
51	Blaydes Lens Holding Forceps very delicate narrow jaws angled 45°	IW-933
52	Dodick Lens Holding Forceps angled smooth lower shaft 12mm & upper shaft 11mm to have a frm grip on iol's	IW-935
53	Kraff Lens Holding Forceps with 3mm curved tips to grasp to the optic & the haptics.	IW-936
54	Kratz Lens Holding Forceps very delicate narrow jaws curved	IW-937
55	Tennant Tying Forceps extra delicate smooth jaws 6mm long tying platform round handle with guide pin	IW-938
56	Mc pherson Tying Forceps 5.00mm long smooth jaws straight with tying platform, big	IW-939
57	Mc pherson Tying Forceps 5.00mm long smooth jaws straight with tying platform, small	IW-939A
58	Mc pherson Tying Forceps curved with smooth jaws, ideal for 8 - 0 to 11 - 0 suture with extra delicated 6mm tying flatform, big	IW-940
59	Mc pherson Tying Forceps curved with smooth jaws, ideal for 8 - 0 to 11 - 0 suture with extra delicated 6mm tying flatform, small	IW-940A

60	Kelmen - Mc pherson Forceps 11.0 mm long platform angeled with smooth jaw ideal for implanting iol & 8 - 0 to 11 - 0 suture, big	IW-941
61	Kelmen - Mc pherson Forceps 11.0 mm long platform angeled with smooth jaw ideal for implanting iol & 8 - 0 to 11 - 0 suture, small	IW-941A
62	Jaffe Tying Forceps with smooth jaws ideal for 8-0 to 11 - 0 sutures due to 6.00 mm long platform	IW-942
63	Mc pherson with smooth jaws 11.00 Mm long platform angled, ideal for implanting IOLs	IW-943
64	Jaffe Capsulorhexis Forceps very delicate round grasping tips extremely thin 11.00mm long shanks	IW-944
65	Utrata Capsulorhexis Forceps very delicate triangular grasping tips extremely thin 11.00mm long shanks	IW-946
66	Utrata Capsulorhexis Forceps very delicate triangular grasping tips extremely thin 11.00mm long shanks round handle with guide pin	IW-947
67	Nevyas Capsulorhexis Forceps cystotomy shaped tips extremely thin 11.00mm long shanks for penetrating and grasping the anterior capsule	IW-948
68	Masket Capsulorhexis Forceps very delicate tringular grasping tips extremely thin 11.00 mm long vaulted shanks	IW-949
69	Inamura Capsulorhexis Forceps very delicate grasping tips and extremely thin 11.00 mm long curveed	IW-950
70	Perone lasik flap forceps abrasive blasted inner surface.	IW-953
71	Jewelers Standard Forceps abrasive blasted inner surface	IW-954
72	Jewelers Delicate Forceps and delicate short forceps	IW-955
73	Akahoshi Phaco Pre Chopper with flattened pointed tips angled and gently curved	IW-956
74	Hartman Hemostatic Mosquito Forceps 20 mm serrated jaws with ratchet lock, 3.5, straight	IW-957
75	Hartman Hemostatic Mosquito Forceps 20 mm serrated jaws with ratchet lock, 3.5, curved	IW-957A
76	Halstead Hemostatic Mosquito Forceps 25 mm serrated jaws with ratchet lock, 5.00, straight & curved	IW-958
77	Indo webal Cheatle's Forceps for lifting instriments from sterilizer overall length 8.0 " & 10.0"	IW-959
78	Indo webal Sponge Holding Forceps serrated jaws overall length 8.0 " & 10.0"	IW-960

